17 beta-Estradiol

CE

Only for in-vitro diagnostic use

Instructions for use / Notice d'utilisation / Instrucciones de uso

English	2
Français	9
Español	16
Literature / Bibliographie / Bibliografía.....	26
Packaging materials / Matériels d'emballage / Materiales de embalaje	26
Symbols Key / Explication des symboles / Símbolos	27
Summary of Test Procedure / Résumé de la procedure de test / Resumen de la técnica	28

REF

DNOV003 (96 Determinations)

ENGLISH

1. INTRODUCTION

17 beta-Estradiol (also known as 17 β -Estradiol, oestradiol or E2) is one of the three main estrogens produced in the human body. The other two are Estrone (E1) and Estriol (E3) which are mainly produced in post-menopausal women and during pregnancy respectively. Estradiol is the predominant estrogen synthesised in pre-menopausal women and has a key role in the regulation of fertility in women.

Estradiol is synthesised in the adrenal glands, the ovaries (women) and tests (men). It is produced by the conversion of testosterone to estradiol by the action of the enzyme aromatase. The enzyme aromatase is also responsible for production of estrone by the conversion of androstenedione¹.

Estrogens regulate the development of secondary sex characteristics; estradiol is also involved in other physiological functions in different organs and systems such as skin, blood vessels, bone, muscle, gastrointestinal tract, brain, lung and pancreas² etc. Quantification of estrogens (and androgens) are useful in the assessment and management of different sex-hormone related disorders including hypogonadism, hirsutism, amenorrhea, fertility and ovarian tumours.

Due to the importance of estradiol in the regulation of the menstrual cycle, low levels are associated with disturbance of the cycle and amenorrhea. Measurements can therefore be used to assess ovarian functionality and support diagnosis of ovarian function recovery in women treated with aromatase inhibitors³.

Estradiol levels are known to be low in men where its biological function is not clearly understood. There is some indication that it is involved in several physiological and pathological conditions such as prostate cancer and infertility. Conversely, elevated levels in men, due to high activity of aromatase, inhibits the synthesis of luteinising hormone from the pituitary gland leading to a decrease in testosterone and onset of hypogonadism and gynecomastia.

Estradiol has a key role in bone homeostasis where low levels are associated with loss of bone mass. It has influence on osteoblast and osteoclast activity to support the maintenance of bone homeostasis and is considered a good predictor of bone mass^{1,2,4-6}.

2. INTENDED USE

17 beta-Estradiol is a manual *in vitro* diagnostic device intended for the quantitative determination of 17 β -Estradiol in human serum or plasma.

3. PRINCIPLE OF THE ASSAY

17 beta-Estradiol is a competitive enzyme immunometric assay (ELISA).

Microtiterplates are coated with anti-17 β -Estradiol antibodies (solid-phase). 17 β -Estradiol in the sample competes with added horseradish peroxidase labelled 17 β -Estradiol (enzyme-labelled antigen) for antibody binding. After incubation, the bound/free separation is performed by a simple solid-phase washing. Then the enzyme HRP in the bound-fraction reacts with the Substrate (H₂O₂) and the TMB Substrate and develops a blue colour that changes into yellow when the Stop Solution (H₂SO₄) is added. The colour intensity is inversely proportional to the 17 β -Estradiol concentration in the sample. Absorption at 450 nm is read using an ELISA Microtiterplate.

17 β -Estradiol concentration in the sample is calculated through a calibration curve.

4. MATERIALS

4.1. Reagents supplied

- **Microtiterplate:** 12 breakapart 8-well snap-off strips coated with anti-17 β -Estradiol antibodies, in resealable aluminium foil.
- **Stop Solution:** 1 bottle containing 15 ml sulphuric acid, 0.15 M (avoid any skin contact).
- **Conjugate:** 1 bottle containing 22 ml of horseradish peroxidase labelled 17 β -Estradiol; contains $\leq 0.0015\%$ (v/v) CMIT/ MIT (3:1) and BSA.
- **TMB Substrate Solution:** 1 bottle containing 15 ml 3, 3', 5, 5'-tetramethylbenzidine (H₂O₂-TMB 0.26g/L) (avoid any skin contact); contains ≤ 0.06 (v/v) CMIT/ MIT (3:1).
- **Wash Solution 10x conc.:** 1 bottle containing 50 ml of a 10x concentrated solution of phosphate buffer 0.2 M, pH 7.4; contains ≤ 0.06 (v/v) CMIT/ MIT (3:1).
- **Control:** 1 bottle containing 0.5 ml of a lot-specific control solution. The concentration is indicated on the label of the bottle; contains ≤ 0.06 (v/v) CMIT/ MIT (3:1) and BSA.
- **Standards:** 6 bottles, each containing 0.5 ml; contain ≤ 0.06 (v/v) CMIT/ MIT (3:1) and human serum.

Standard 0:	0 pg/ml
Standard 1:	20 pg/ml
Standard 2:	120 pg/ml
Standard 3:	300 pg/ml
Standard 4:	600 pg/ml
Standard 5:	2000 pg/ml

For hazard and precautionary statements see 16.1.

4.2. Materials supplied

- 1 Cover foil
- 1 Instructions for use (IFU)

4.3. Materials and Equipment needed

- ELISA Microtiterplate reader, equipped for the measurement of absorbance at 450/620-630 nm
- Incubator 37 °C
- Manual or automatic equipment for rinsing Microtiterplates
- Pipettes to deliver volumes between 25 and 1000 µL
- Vortex tube mixer
- Distilled water
- Disposable tubes

5. STABILITY AND STORAGE

The reagents are stable up to the expiry date stated on the label when stored at 2...8 °C in the dark. Once opened, the kit is stable at 2...8°C for 6 months.

6. REAGENT PREPARATION

It is very important to bring all reagents, samples and standards to room temperature (22...28 °C) before starting the test run for at least 30 minutes! At the end of the assay, store immediately the reagents at 2...8 °C; avoid long exposure to room temperature.

6.1. Microtiterplate

The break-apart snap-off strips are coated with anti-17β-Estradiol antibodies. Store at 2...8 °C. Open the bag only when it is at room temperature. Immediately after removal of strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C; stability until expiry date.

6.2. Conjugate

The conjugate is ready to use. Mix gently on a rotating mixer for 5 min.

6.3. Standards

The standard solutions are ready to use.

6.4. TMB Substrate Solution

The bottle contains 15 ml of a tetramethylbenzidine/hydrogen peroxide system. The reagent is ready to use and has to be stored at 2...8 °C in the dark. The solution should be colourless or could have a slight blue tinge. If the substrate turns into blue, it may have become contaminated and should be thrown away.

6.5. Stop Solution

The bottle contains 15 ml 0.15 M sulphuric acid solution (H314). This ready to use solution has to be stored at 2...8 °C.

6.6. Wash Solution

Dilute the concentrated solution with distilled water to a final volume of 500 ml prior to use. For smaller volumes respect the 1:10 dilution ratio. The diluted solution is stable for 30 days at 2...8 °C. In the concentrated solution it is possible to observe the presence of crystals, in this case mix at room temperature until complete dissolution of crystals. For greater accuracy dilute the whole bottle of concentrated wash solution to 500 ml and take care that all crystals are transferred by washing the bottle, then mix until crystals are completely dissolved.

6.7. Control

The bottle contains 0.5 ml of a lot-specific control solution. The concentration is indicated on the label. The control is ready to use.

7. SPECIMEN COLLECTION AND PREPARATION

The assay should be performed using serum (standard sampling tubes or tubes containing serum separating gel) or plasma (lithium heparin, sodium heparin or potassium EDTA) samples. Store the sample at -20°C if the determination is not performed on the same day of the sample collection. Avoid repetitive freezing and thawing of samples. Before using, mix gently for 5 minutes with a rotating mixer.

<i>Sample Storage</i>	<i>Duration</i>
2...8 °C	24 hours
Freeze/thaw cycles	1 cycle

8. ASSAY PROCEDURE

8.1. Test Preparation

Please read the Instructions for use carefully **before** performing the assay. Result reliability depends on strict adherence to the Instructions for use as described. Prior to commencing the assay, the distribution and identification plan for all specimens and standards should be carefully established. Select the required number of microtiter strips or wells and insert them into the holder.

Please allocate at least:

1 well	(e.g. A1)	for the substrate blank
2 wells	(e.g. B1+C1)	for standard 0
2 wells	(e.g. D1+E1)	for standard 1
2 wells	(e.g. F1+G1)	for standard 2
2 wells	(e.g. H1+A2)	for standard 3
2 wells	(e.g. B2+C2)	for standard 4
2 wells	(e.g. D2+E2)	for standard 5
2 wells	(e.g. F2+G2)	for control

It is necessary to determine standards, control and patient samples in duplicate.

Perform all assay steps in the order given and without any appreciable delays between the steps.

A clean, disposable tip should be used for dispensing each standard and each patient sample.

Adjust the incubator to 37 °C.

1. Dispense 25 µl standards, control and samples into their respective wells. Add 200 µl conjugate to each well. Leave well A1 for substrate blank.
2. Cover wells with the foil supplied in the kit.
3. **Incubate for 2 hour at 37 °C.**
4. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 µl diluted washing solution. Avoid overflows from the reaction wells.

Important note: During each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.

Automatic washer: In case you use automatic equipment, wash the wells at least 5 times.

Note: Washing is critical! Insufficient washing results in poor precision and falsely elevated absorbance values.

5. Dispense 100 µl TMB Substrate Solution into all wells.
6. **Incubate for exactly 30 min at room temperature (22...28 °C) in the dark.**
7. Dispense 100 µl Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate Solution. Shake the microplate gently.
Any blue colour developed during the incubation turns into yellow.
8. Measure the absorbance (E) of the specimen at 450 nm against a reference wavelength of 620-630 nm or against blank within 5 minutes.

8.2. Measurement

Adjust the ELISA Microwell Plate Reader **to zero** using the **substrate blank in well A1**.

If - due to technical reasons - the ELISA reader cannot be adjusted to zero using the substrate blank in well A1, subtract the absorbance value of well A1 from all other absorbance values measured in order to obtain reliable results!

Measure the absorbance of all wells at **450 nm** and record the absorbance values for each standard and patient sample.

Where applicable calculate the **mean absorbance values** of all duplicates.

9. QUALITY CONTROL

Good Laboratory Practice (GLP) requires the use of quality control specimens in each series of assays in order to check the performance of the assay. Controls should be treated as unknown samples, and the results analysed with appropriate statistical methods.

The kit control provided in the kit should be tested as unknown and is intended to assist in assessing the validity of results obtained with each assay plate.

We recommend the users to maintain graphic records of the control values generated with each assay run, including the running means, SDs and %CVs. This information will facilitate the controls trending analysis relating to the performance of current and historical control lots relative to the supplied Quality Control data. The trending will assist in the identification of assays which give control values significantly different from their average range.

When interpreting control data, users should note that this product was designed and developed as a manual product. The range stated on the Certificate of Analysis (CoA) should be appropriate for assays that are performed manually and with strict adherence to the Assay Procedure described above. It is recognised by Quality Control professionals, that as a result of differences in conditions and practices, there will always be variability in the mean values and precision of control measurements between different laboratories⁷.

10. CALCULATION OF RESULTS

A variety of data reduction software packages are available, which may be employed to generate the mean calibration curve and to calculate the mean concentrations of unknown samples and controls. A 4-parameter logistic (4PL) curve fit, **including Standard 0 is recommended.**

Alternatively, a calibration curve may be prepared on semi-log graph paper by plotting mean absorbance on the Y-axis against concentration of analyte on the X-axis. Standard 0 should be included in the calibration curve. Read the mean absorbance value of each unknown sample off the curve.

Conversion of units

To convert results to SI units:

pmol/L = pg/mL x 3.671

11. MEASURING RANGE

The assay measuring range (AMR) is 29.6 – 2000 pg/mL (108.7 – 7342 pmol/L).

Any value that reads below 29.6 pg/mL (108.7 pmol/L) should be reported as “< 29.6 pg/mL (108.7 pmol/L)”. Any value that reads above 2000 pg/mL (7342 pmol/L) should be reported as “> 2000 pg/mL (7342 pmol/L)”.

12. METROLOGY AND TRACEABILITY

The standards of this kit are traceable to the Estradiol standard from the National Metrology Institute of Japan (6004-a).

13. EXPECTED VALUES

The following ranges were determined using the 17 beta-Estradiol and are provided for information only. The 90 % reference interval for apparently healthy adults were calculated by a non-parametric method following guidance from CLSI C28-A3 “Defining, Establishing and Verifying Reference Intervals in the Clinical Laboratory”.

Adults	n	Median (pg/mL)	Reference interval (pg/mL)
Males	120	34.3	<29.6 – 64.2
Females			
Luteal phase	40	71.6	<29.6 – 212.4
Follicular phase	39	46.5	<29.6 – 135.6
Ovulatory phase	41	57.1	<29.6 – 273.5
Post-menopausal	120	32.9	<29.6 – 65.8
Pregnancy			
1 st trimester	40	1120.1	301.2 – >2000
2 nd trimester	39	>2000	1850.3 – >2000
3 rd trimester	40	>2000	>2000

Children	No. of subjects	Median (pg/mL)	Reference interval (pg/mL)
Females			
1 – 10 years	39	56.7	<29.6 – 115.1
10 – 15 years	41	79.8	46.9 – 155.5
15 – 19 years	40	138.3	88.8 – 264.1
Males			
1 – 10 years	40	58.5	<29.6 – 115.7
10 – 15 years	39	39.2	<29.6 – 105.1
15 – 19 years	39	31.7	<29.6 – 91.3

The above ranges should be considered as guidelines only; it is recommended that each laboratory establish its own expected range based upon its own patient population.

14. SPECIFIC PERFORMANCE CHARACTERISTICS

Representative performance data are shown. Results obtained at individual laboratories may vary.

14.1. Detection Capability

The limit of blank (LoB), limit of detection (LoD) and limit of quantitation (LoQ) were determined with guidance from CLSI EP17-A, “Protocols for Determination of Limits of Detection and Limits of Quantitation” using 6 blanks and 6 low level samples.

Sensitivity	Concentration
Limit of Blank (LoB)	6.8 pg/mL
Limit of Detection (LoD)	14.6 pg/mL
Limit of Quantitation (LoQ)	29.6 pg/mL

14.2. Trueness

Trueness has been demonstrated through a recovery test for the 17 beta-Estradiol with the 17β-Estradiol standard from the National Metrology Institute of Japan (code 6004-a).

14.3. Precision

Precision of the 17 beta-Estradiol was determined by performing a complex precision study.

Repeatability

A total of 6 serum samples were assayed in 5 replicates, once a day for 5 days by 3 operators. Data from one representative lot is shown below:

Sample	n	Mean Conc. (pg/mL)	Within run (Repeatability)	
			SD	CV%
1	75	51.7	6.1	11.9%
2	75	134.2	11.3	8.4%
3	75	264.0	25.4	9.6%
4	75	627.7	56.4	9.0%
5	75	838.1	61.5	7.3%
6	75	1346.6	120.2	8.9%

Reproducibility

A total of 6 serum samples were assayed in 5 replicates, once a day for 5 days by 3 operators. Results for the combined data from 2 lots is shown below:

Sample	n	Mean Conc. (pg/mL)	Within Laboratory (Reproducibility)	
			SD	CV%
1	150	52.4	8.2	15.7%
2	150	139.0	21.2	15.3%
3	150	274.7	39.0	14.2%
4	150	633.9	66.0	10.4%
5	150	860.0	95.3	11.1%
6	150	1374.7	159.8	11.6%

14.4. Linearity

Linearity was evaluated based on CLSI EP-06, "Evaluation of the Linearity of Quantitative Measurement Procedures". For estradiol concentration by 17 beta-Estradiol, the measurement procedure shows linearity for the interval from 16.5 to 2271.3pg/mL (60.57 to 8337.9 pmol/L) within the allowable deviation of linearity (ADL) of $\pm 15\%$.

14.5. Method comparison

The 17 beta-Estradiol was compared against a commercially available quantitative manual assay, following CLSI EP-09C, "Method Comparison and Bias Estimation Using Patient Samples". A total of 106 samples, selected to represent a wide range of estradiol concentrations, was assayed by each method.

Passing-Bablok regression analysis was performed on the comparative data:

n	Slope [95% CI]	Intercept (pg/mL) [95% CI]	Correlation coefficient (r)
106	0.98 [0.93 to 1.04]	14.5 [8.6 to 18.2]	0.99

14.6. Analytical Specificity

The specificity was assessed with the following cross-reactants.

Cross reactant	Concentration tested	Mean % Cross reactivity
Estriol	100 ng/mL	1.0%
Estrone	100 ng/mL	0.9%
Fulvestrant	100 ng/mL	0.4%
Testosterone	100 ng/mL	0.0%
Cortisol	1000 ng/mL	0.0%
Progesterone	200 ng/mL	0.0%
DHEA-S	10000 ng/mL	0.0%
17OH Progesterone	100 ng/mL	0.0%
5 α -dihydrotestosterone (DHT)	500 ng/mL	0.0%
Androstenedione	100 ng/mL	0.0%
Androsterone	100 ng/mL	0.0%

Cross reactant	Concentration tested	Mean % Cross reactivity
Cortisone	1000 ng/mL	0.0%
Prednisolone	200 ng/mL	0.0%
Danazol	1000 ng/mL	0.0%
2-Methoxyestradiol	10 ng/mL	0.7%
17 α -Ethinylestradiol	200 ng/mL	0.0%
17 β -Estradiol 17-sulphate	100 ng/mL	0.0%

The following substances do not interfere with a bias of $> \pm 15\%$ in the 17 beta-Estradiol assay when the concentrations are below the stated threshold presented in the following table.

Potentially Interfering Reagent	Threshold Concentration
Bilirubin, conjugated	3.89 mg/dL
Bilirubin, unconjugated	20 mg/dL
Haemoglobin	640 mg/dL
Total Protein	8.3 g/dL
Triglycerides	2000 mg/dL

14.7. Serum-plasma study

The 17 beta-Estradiol matrix comparison study was performed to evaluate the difference across tube types (serum separator tubes (SST), lithium heparin plasma, sodium heparin plasma and K2 EDTA plasma) versus the control samples (red top serum, without additive) following CLSI EP-35 guidelines. A total of 22 samples (18 native, 4 spiked) to cover the range were evaluated. Passing-Bablok regression analysis was performed on the comparative data:

Sample type	Slope [95% CI]	Intercept (pg/mL) [95% CI]	Correlation coefficient (r)
SST	1.03 [0.98 to 1.07]	-3.5 [-37.2 to 30.2]	1.00
Lithium Heparin	1.04 [0.97 to 1.12]	-0.8 [-59.0 to 57.4]	0.99
Sodium Heparin	1.10 [1.03 to 1.17]	-8.05 [-62.2 to 46.1]	0.99
EDTA	1.06 [1.01 to 1.11]	-7.6 [-46.3 to 31.1]	0.99

15. LIMITATIONS OF USE

- As in the case of any diagnostic procedure, results must be interpreted in conjunction with the patient's clinical presentation and other information available to the physician.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays⁸. Patients routinely exposed to animals or to animal serum products can be prone to this interference and anomalous values may be observed.

16. PRECAUTIONS AND WARNINGS

- The test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The use of the testkits with analyzers and similar equipment has to be validated. Any change in design, composition and test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents for these reasons. The manufacturer is not liable for any results by visual analysis of the patient samples.
- This kit is intended for *in vitro* use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- All human source material used in the preparation of the reagents has been tested and found negative for antibody to HIV 1&2, HbsAg, and HCV. No test method however can offer complete assurance that HIV, HBV, HCV or other infectious agents are absent. Therefore, the Calibrators and the Controls should be handled in the same manner as potentially infectious material.
- Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.
- Some reagents contain small amounts of Sodium Azide (NaN₃) as preservative.
- Sodium Azide may be toxic if ingested or absorbed through the skin or eyes; moreover, it may react with lead or copper plumbing to form potentially explosive metal azides. If you use a sink to remove the reagents, wash through large with amounts of water to prevent azide build-up.
- The TMB Substrate contains an irritant, which harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.

- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous, corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H₂O₂ to direct sunlight, metals or oxidants. Do not freeze the solution.
- The clinical significance of the determination 17β-Estradiol can be invalidated if the patient was treated with cortisone or natural or synthetic steroids.
- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this instructions for use.
- All reagents should be stored refrigerated at 2...8°C in their original container. Any exceptions are clearly indicated.
- Allow all kit components and specimens to reach room temperature (22...28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate.
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.

16.1. Safety note for reagents containing hazardous substances

Reagents may contain CMIT/MIT (3:1) (refer to 4.1).

Therefore, the following hazard and precautionary statements apply.

	Warning	
	H317	May cause an allergic skin reaction.
	P261	Avoid breathing dust/fumes/gas/mist/vapours/spray.
	P280	Wear protective gloves/protective clothing/eye protection/face protection.
	P302+P352	IF ON SKIN: Wash with plenty of water.
	P333+P313	If skin irritation or rash occurs: Get medical advice/attention.
	P362+P364	Take off contaminated clothing and wash it before reuse.

Further information can be found in the safety data sheet.

16.2. Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

For information about the packaging materials refer to PACKAGING MATERIALS.

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
LITERATURE / BIBLIOGRAPHIE / BIBLIOGRAFÍA

1. Smy L, Straseski JA. Measuring estrogens in women, men, and children: Recent advances 2012-2017. *Clin Biochem.* 2018 Dec;62:11-23.
2. Rosner W, Hankinson SE, Sluss PM, Vesper HW, Wierman ME. Challenges to the measurement of estradiol: an endocrine society position statement. *J Clin Endocrinol Metab.* 2013 Apr;98(4):1376-87.
3. van Hellemond IEG, Vriens IJH, Peer PGM, Swinkels ACP, *et al.*; Dutch Breast Cancer Research Group (BOOG). Ovarian Function Recovery During Anastrozole in Breast Cancer Patients With Chemotherapy-Induced Ovarian Function Failure. *J Natl Cancer Inst.* 2017 Dec 1;109(12).
4. Katulski K, Slawek S, Czyzyk A, Podfigurna-Stopa A, *et al.* Bone mineral density in women with polycystic ovary syndrome. *J Endocrinol Invest.* 2014 Dec;37(12):1219-24.
5. Cameron DA, Douglas S, Brown JE, Anderson RA. Bone mineral density loss during adjuvant chemotherapy in premenopausal women with early breast cancer: is it dependent on oestrogen deficiency? *Breast Cancer Res Treat.* 2010 Oct;123(3):805-14.
6. Decaroli MC, Rochira V. Aging and sex hormones in males. *Virulence.* 2017 Jul 4;8(5):545-570.
7. Basic QC Practices On-line Course; <http://www.Westgard.com>.
8. Boscato, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27–33

PACKAGING MATERIALS / MATERIELS D'EMBALLAGE / MATERIALES DE EMBALAJE

 PAP 21	 PAP 21	 PAP 22	MTP  ALU / LDPE 90
---	---	--	--

SYMBOLS KEY / EXPLICATION DES SYMBOLES / SIMBOLOS

	Manufactured by / Fabriqué par/ Fabricado por
IVD	In Vitro Diagnostic Medical Device / Dispositif médical de diagnostic <i>in vitro</i> / Producto para diagnóstico In vitro
LOT	Lot Number / Numéro de lot / Número de lote
	Expiration Date / Date de péremption / Fecha de caducidad
	Storage Temperature / Température de conservation / Temperatura de almacenamiento
	Keep away from sunlight / Protéger de rayonnement solaire / Mantener alejado de la luz solar
CE	CE Mark / Marquage CE / Marca CE
REF	Catalogue Number / Référence du catalogue / Número de Catálogo
	Consult Instructions for Use / Consulter la notice d'utilisation / Consulte las Instrucciones de Uso
MTP	Microtiterplate / Plaque de Microtitrage / Placa de Microtitulación
CONJ	Conjugate / Conjugué / Conjugado
CONTROL	Control / Contrôle / Control
CAL	Calibrator resp. Standard / Standard o Étalon / Estándar o Calibrador
SOLN STOP	Stop Solution / Solution d'arrêt / Solución de Parada
SUB TMB	TMB Substrate Solution / Solution de Substrat TMB / Solución de Sustrato TMB
WASH BUF 10x	Wash Solution 10x concentrated / Solution de lavage concentré 10 x / Solución de lavado concentrado x10
	Contains sufficient for "n" tests/ Contenu suffisant pour "n" tests / Contenido suficiente para "n" tests

SCHEME OF THE ASSAY

17 beta-Estradiol

Test Preparation

Prepare reagents and samples as described.
Establish the distribution and identification plan for all samples and standards/controls.
Select the required number of microtiter strips or wells and insert them into the holder.

Assay Procedure

	Substrate blank	Standard 1 – 5	Control	Sample
Standard 1 – 5	-	25 µl	-	-
Control	-	-	25 µl	-
Sample	-	-	-	25 µl
Conjugate	-	200 µl	200 µl	200 µl
Cover wells with foil supplied in the kit Incubate for 2 hour at 37 °C Wash each well three times with 300 µl diluted wash solution In case you use automatic equipment, wash the wells at least 5 times.				
TMB Substrate Solution	100 µl	100 µl	100 µl	100 µl
Incubate for exactly 30 min at room temperature (22...28 °C) in the dark				
Stop Solution	100 µl	100 µl	100 µl	100 µl
Shake the microplate gently Photometric measurement at 450 nm, 620 – 630 nm				

Gold Standard Diagnostics Frankfurt GmbH

Waldstrasse 23 A6

63128 Dietzenbach, Germany

Tel.: +49 6074 23698-0

Fax: +49 6074 23698-900

E-Mail: info.frankfurt@eu.goldstandarddiagnostics.com

Website: clinical.goldstandarddiagnostics.com

DNOV003_Estradiol_IFU_rev02_fromLot_6090BN