GSD NovaPrime[®] TSP SARS-CoV-2

RT-PCR

Only for in-vitro diagnostic use

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REF

PCOV6111
PCOV6113

(1 x 96 Determinations)
(3 x 96 Determinations)

ENGLISH

1. INTENDED USE

The GSD NovaPrime® TSP SARS-CoV-2 RT-PCR is intended for the qualitative determination of SARS-CoV-2 (severe acute respiratory syndrome coronavirus 2) genomic RNA extracted from human upper respiratory (nasal wash/swab, nasopharyngeal swab, oropharyngeal swab, saliva, and gargling solution) specimen types collected from individuals suspected of respiratory viral infection.

2. PRINCIPLE OF THE ASSAY

The qualitative determination of specific RNA is based on Real-Time reverse-transcription Polymerase Chain Reaction (RT-PCR) technology. The kit contains specific primers and probes labelled with fluorescent reporter and quencher dyes for amplification and simultaneous detection of RNA sequences which represent two specific regions of the SARS-CoV-2 N gene. Furthermore, the assay contains a heterologous amplification target (Extraction Control, **EC**) to identify possible RT-PCR inhibition by interfering substances contained in the sample or failure of the preceding RNA extraction. Therefore, the **EC** is added to the specimen during RNA isolation.

The gene of interest specific probes are labelled with the fluorophore FAM™. The **EC** specific probe is labelled with the fluorophore Cy5™ thereby allowing parallel detection of both amplicons in the corresponding detection channels.

The GSD NovaPrime® TSP SARS-CoV-2 has been validated for following Real-Time PCR platforms:

- AriaDx™ (Agilent Technologies)
- AriaMx™ (Agilent Technologies)
- ABI Prism® 7500 Fast SDS (Applied Biosystems™) in Standard and Fast Mode
- CFX96 Touch™ (Bio-Rad)
- LightCycler® 480 II (Roche)
- QuantStudio™ 5 (Applied Biosystems™)

3. MATERIALS

3.1. Reagents Provided

Cap	Symbol	Component	Volume per Vial [µL]	Number of Vials PCOV6111 (PCOV6113)
green	POL 2x	Polymerase (hot-start DNA polymerase, nucleotides, magnesium, enhancers and stabilizers)	1060	1 (3)
green/white	RT 20x	reverse transcriptase with RNase inhibitor	106	1 (3)
blue	PP TSP	Primer-Probe-Mix	106	1 (3)
red	PC	Positive Control (plasmid DNA representing the N gene of SARS-CoV-2)	150	1 (3)
yellow	EC	Extraction Control (RNA-based internal lysis, extraction, and amplification control)	500	1 (3)
transparent	NFW	Nuclease Free Water	500	1 (3)

3.2. Materials and Equipment needed, but not provided

- Biological safety cabinet for sample handling
- Biohazard waste container
- Equipment and consumables for isolating virus RNA from respiratory specimens:
Appropriate nucleic acid extraction kits, e.g., Mag-Bind® Viral DNA/RNA Xpress Kit (Omega Bio-tek Inc., Cat. No. M6219-2304CEIVD) validated for Hangzhou Allsheng Instruments Co. Ltd. Auto-Pure96.
Use of other RNA extraction methods must be validated by the user.
- Real-Time PCR instrument (for already validated instruments refer to 2. PRINCIPLE OF THE ASSAY). Alternative Real-Time PCR instruments might also be appropriate. Their suitability for use with the GSD NovaPrime® TSP SARS-CoV-2 has to be validated by the user.
- Appropriate Real-Time PCR consumables (e.g., disposable tubes, reaction plates, corresponding optical closing materials)
- Benchtop microcentrifuge
- Centrifuge with a rotor for microtiter plates
- Vortex mixer
- Adjustable pipettes in relation to reaction setup
- Disposable DNase/RNase free pipette tips with filters
- Disposable powder-free gloves

4. STABILITY AND STORAGE

The GSD NovaPrime® TSP SARS-CoV-2 kit is shipped on dry ice and all components should arrive frozen.

- All components have to be stored between -30 °C and -15 °C immediately after arrival.
- Repeated freeze thaw cycles (more than seven) of reagents should be avoided since this might affect the performance of the kit. Reagents should be frozen in aliquots if they are used intermittently.
- Keep unfrozen storage (e.g., storage on ice) as short as possible.
- [POL 2x], [RT 20x] and [PP TSP] should be stored frozen and protected from light until use.

5. SAMPLE PREPARATION

- Extracted RNA or total nucleic acid extracted from human upper respiratory (nasal wash/swab nasopharyngeal swab, oropharyngeal swab, saliva, and gargling solution) specimen types is the starting material for the GSD NovaPrime® TSP SARS-CoV-2.
- The quality and the quantity of the extracted RNA has a crucial effect on the performance of the entire RT-PCR test system. Make sure that the nucleic acid extraction method is compatible with Real-Time PCR technology.
- Make sure that sampling conditions (time, method, sample storage and transport) are suitable, to avoid negative impact on the extraction and RT-PCR.
- For nucleic acid extraction a method suitable for extracting virus RNA from human respiratory specimen should be used.
- The [EC] can be used for monitoring both the RNA extraction procedure and any potential PCR inhibition. Therefore, the [EC] has to be added prior to the nucleic acid extraction procedure into the lysis buffer. Independent of the method / system used for nucleic acid extraction, the [EC] can be directly added to the specimen.
- About 5 µL [EC] may be suitable but should be carefully evaluated. Excess amounts of [EC] may lead to a reduced amplification of the SARS-CoV-2 specific sequences and might therefore result in increased FAM™ Ct values. Underestimation of the SARS-CoV-2 specific signal may occur.
- Since ethanol is a strong Real-Time PCR inhibitor, it is necessary to completely eliminate it prior to the elution of the nucleic acid during extraction. If using spin columns with washing buffers containing ethanol, it is highly recommended to perform an additional centrifugation step of 10 min at approximately 17000 x g (~ 13000 rpm) before eluting the RNA. For this additional centrifugation step, use a new collection tube.
Concentrations > 4 % ethanol (v/v) resulted in inhibition of the GSD NovaPrime® TSP SARS-CoV-2 assay.

6. ASSAY PROCEDURE

6.1. Reaction Setup

- Please read the instructions for use carefully before performing the assay. Reliability of results depends on following strictly the instructions for use.
- Before use make sure that all samples and reagents are thawed completely, mixed by up and down pipetting or vortexing and centrifuged briefly.
- **Important:** [RT 20x] and [POL 2x] are viscous. Briefly spin down the tubes to ensure material has not lodged in the cap or side of tube. Be sure to pipette and dispense carefully and use pipette tips suitable for pipetting viscous liquids. Keep all components ([POL 2x], [RT 20x] and [PP TSP]) permanently on ice/cooling block during PCR setup. Prepare final reaction mix fresh each time and immediately before starting the Real-Time RT-PCR run.
- The use of [NFW] as no template control (NTC) is highly recommended.
- Define the positions (wells) for samples and controls ([PC], NTC) on the plate.

Reaction Setup	
Component	Volume
[POL 2x]	10 µL
[RT 20x]	1 µL
[PP TSP]	1 µL
Sample or [PC] or [NFW]	8 µL
Total volume	20 µL

- Close the optical reaction plate with corresponding optical closing material.
- Centrifuge the optical reaction plate in a centrifuge with a rotor for microtiter plates for 60 seconds at 4 °C at approximately 1000 x g (~ 3000 rpm).

6.2. Programming the Real-Time PCR Instrument

Regarding setup and programming of the Real-Time PCR instrument, please use the manual of the respective instrument.

For detailed programming instructions regarding the use of the GSD NovaPrime® TSP SARS-CoV-2 on specific Real-Time PCR instruments please contact the manufacturer.

RT-PCR Run Settings	
Reaction Volume	20 µL
Passive Reference	ROX™

Detection Channels

Detection of the amplified viral nucleic acid fragments is performed in the following detection channels:

Target	Fluorophore (Quencher)	Detection Channel				
		AriaDx™ AriaMx™	CFX96 Touch™	LightCycler® 480 II	ABI Prism® 7500 SDS/Fast SDS	QuantStudio™ 5
SARS-CoV-2	FAM™ (BHQ1)	FAM	FAM	FAM (465-510)	FAM	x1-m1
EC	Cy5™ (BHQ2)	Cy5	Cy5	Cy5 (618-660)	Cy5	x5-m5

Temperature Profile and Data Collection

No. of Cycles	Temperature	Time (min)	Data Collection
1	45 °C	10:00	-
1	95 °C	02:00	-
40*	95 °C	00:05	-
	60 °C	00:30	Fluorescence measurement at the end of every cycle

* LC480 II: 45 cycles

Before starting the test run, please check the settings for cycles, temperature, and time.

7. RESULTS

Data analysis should be performed with the software of the used real-time PCR device according to the manufacturer's instructions.

Analysis settings

RT-PCR instrument	Analysis	Baseline
AriaDx™, AriaMx™ (Agilent Technologies)	Threshold 0.1	Select start and end cycle values so that possible initial signal noise is not included, and the baseline ends before significant fluorescence is detected.
ABI Prism® 7500 SDS Fast (Applied Biosystems™)		
QuantStudio™ 5 (Applied Biosystems™)		
CFX96 Touch™ (Bio-Rad)	Auto Threshold	
LightCycler® 480 II (Roche)	Abs Quant/2 nd Derivative Max	

7.1. Interpretation of Results

The test run is valid only if RT-PCR run is complete.

The GSD NovaPrime® TSP SARS-CoV-2 test protocol dictates that the controls must be analyzed before patient sample results.

7.1.1. Controls

If kit control(s) fail, the test is invalid and needs to be repeated.

Sample data is analyzed and interpreted only after all the kit controls pass.

- NTC (no template control): must NOT have a detectable Ct for any target. If this control has a detectable Ct, this indicates contamination of the PCR run. It is considered invalid and must be repeated.
- PC**: FAM Ct ≤ 38
- EC**: all human samples should exhibit Cy5 Ct values ≤ 35

Failure to detect **EC** in any specimen may indicate:

- Improper extraction of nucleic acid from clinical materials resulting in loss of RNA and/or RNA degradation.
- (RT-) PCR inhibition.
- Improper assay set up and execution.
- Reagent or equipment malfunction.

If the **EC** assay does not produce a positive result for specimens, interpret as follows:

- If the SARS-CoV-2 target is positive even in the absence of a positive **EC** signal, the result should be considered valid. A negative **EC** signal does not preclude the presence of SARS-CoV-2 RNA in a specimen.
- If the SARS-CoV-2 target AND **EC** are negative for the specimen, the result should be considered invalid for the specimen. If residual specimen is available, repeat the extraction procedure and repeat the test. If all markers remain negative after re-test, report the results as invalid and a new specimen should be collected if possible.

7.1.2. Samples

If test run is valid interpretation of sample results is as follows:

- Positive (POS): Ct ≤ 38
- Negative (neg): Not detected (ND) or Ct > 38

Please refer to the following table for interpretation of signal combinations.

Detection Channel		Interpretation of Results
FAM (SARS-CoV-2 N gene)	Cy5 EC	
POS	POS	The sample contains SARS-CoV-2 specific RNA.
POS	neg	The sample contains SARS-CoV-2 specific RNA.
neg	POS	The sample does not contain detectable amounts of SARS-CoV-2 specific RNA.
neg	neg	PCR inhibition or reagent failure. A diagnostic statement must not be made. The RT-PCR run should be repeated, or a new sample must be analyzed.

8. SPECIFIC PERFORMANCE CHARACTERISTICS

The results refer to the groups of samples investigated; these are not guaranteed specifications.

8.1. Clinical Agreement

The determinations were performed on the AriaDx™, AriaMx™ (Agilent Technologies), ABI Prism® 7500 SDS Fast and QuantStudio™ 5 (Applied Biosystems™) Real-Time PCR instruments.

A total of 600 nasopharyngeal swab samples were tested. This comprised:

- 100 samples positive for SARS-CoV-2 RNA (including different variants)
- 500 samples negative for SARS-CoV-2

The results obtained with the GSD NovaPrime® TSP SARS-CoV-2 were compared with the results of a CE marked reference RT-PCR assay.

		Reference (RT-PCR)	
		positive	negative
GSD NovaPrime® TSP SARS-CoV-2	positive	100	0
	negative	0	500

Positive Percent Agreement: 100 %

Negative Percent Agreement: 100 %

Overall Percent Agreement: 100 %

8.2. Precision

Repeatability was analyzed with three samples of different target analyte RNA concentrations (high: 10,000 cp/μL, medium: 500 cp/μL, and low: 3x LoD) measured in 5 replicates in two independent runs per day for 20 days. Reproducibility was determined by analyzing the three samples in 5 replicates by three operators at three sites for 5 days and one run per day. The inter-lot reproducibility was analyzed using the three samples with three GSD NovaPrime® TSP SARS-CoV-2 batches in two runs and nine replicates per lot.

Sample Concentration	SARS-CoV-2 (FAM)								
	Repeatability			Reproducibility			Inter-Lot Reproducibility		
	Mean [Ct]	SD [Ct]	% CV	Mean [Ct]	SD [Ct]	% CV	Mean [Ct]	SD [Ct]	% CV
High	20.97	0.54	2.6	21.12	0.97	4.6	22.09	0.33	1.5
Medium	24.48	0.43	1.8	25.31	1.24	4.9	26.17	0.47	1.8
Low	34.07	0.97	2.8	35.78	4.04	5.7	35.57	1.58	4.4

8.3. Analytical Sensitivity

The Limit of Detection (LoD) is defined as the lowest concentration of target analyte that can be detected with a defined level of confidence. The determinations were performed on the AriaDx™ (Agilent Technologies) Real-Time PCR instrument.

LoD of in-vitro transcribed RNA (IVT-RNA)

Pooled negative nasopharyngeal swab material was extracted and spiked with quantified SARS-CoV-2 N gene at different concentrations. The provisional LoD was established by the lowest concentration that gives 100 % positive results when tested in triplicates. The final LoD was determined using a dilution series around the provisional LoD (Probit Analysis) and was tested in three independent runs with 8 replicates each on 3 lots.

	LoD based on Probit analysis
AriaDx™ (Agilent Technologies)	1.75* copies/reaction

* Based on the Poisson distribution it should be assumed that less than 3 copies of target per PCR reaction cannot be reliably detected. Concentration value refers to the sample after RNA extraction.

LoD of WHO International Standard or inactivated full virus particles

Negative nasopharyngeal matrix was spiked with a dilution series (triplicates) of either First WHO International Standard for SARS-CoV-2 RNA (NIBSC code: 20/146) or chemically inactivated intact SARS-CoV-2 virus particles (NATtrol™ SARS-CoV-2). Extraction was performed using the GSD NovaPrime® IVD RNA Extraction AE1 kit. The provisional LoD was established by the lowest concentration that gives 100 % positive results when tested in triplicates. The final LoD was determined using a dilution series around the provisional LoD (Probit Analysis) and was tested in three independent runs with 8 replicates each.

Analyte	LoD based on Probit analysis**
WHO Standard (NIBSC code: 20/146)	52.33 IU/mL
NATtrol™ SARS-CoV-2	37.46 copies/mL

** Concentration values refer to the sample before RNA extraction

LoD on additional RT-PCR instruments

The LoD for additional RT-PCR instruments was performed by preparing three dilutions around the LoD. Therefore, nasopharyngeal swab material was extracted, spiked with quantified SARS-CoV-2 N gene IVT-RNA, and tested on the respective RT-PCR instrument with 24 replicates. LoD was confirmed (within 2-3-fold) on AriaMx™ (Agilent Technologies), ABI Prism® 7500 SDS Fast (Applied Biosystems™) in Standard and Fast Mode, QuantStudio™ 5 (Applied Biosystems™), CFX96 Touch™ (BioRad) and LightCycler® 480 II (Roche).

8.4. Inclusivity

The analytical sensitivity of the GSD NovaPrime® TSP SARS-CoV-2 is first and foremost ensured by the thorough selection of the oligonucleotides. Inclusivity was evaluated by *in-silico* analysis using publicly available complete SARS-CoV-2 virus sequences from the Global Initiative on Sharing Avian Influenza Data (GISAID) database on October 12th, 2023. Based on the current analysis, all of the analyzed SARS-CoV-2 virus variants are detectable by the GSD NovaPrime® TSP SARS-CoV-2.

The occurrence of single isolates that are not detected by the GSD NovaPrime® TSP SARS-CoV-2 cannot be ruled out due to random mutations in the viral genomes. However, only if a mutation is dominant and replaces other virus variants, this is judged as critical for the performance of the assay.

8.5. Cross-Reactivity

Cross-reactivity was evaluated by *in-silico* analysis against normal flora or pathogens that cause similar symptoms or pathogens related to SARS-CoV-2 virus.

Human coronavirus 229E	Influenza A virus	<i>Streptococcus pneumoniae</i>
Human coronavirus OC43	Influenza B virus	<i>Streptococcus pyogenes</i>
Human coronavirus HKU1	Enterovirus	<i>Bordetella pertussis</i>
Human coronavirus NL63	Respiratory syncytial virus	<i>Mycoplasma pneumoniae</i>
SARS coronavirus	Rhinovirus	<i>Pneumocystis jirovecii</i>
MERS coronavirus	<i>Chlamydia pneumoniae</i>	<i>Candida albicans</i>
Adenovirus	<i>Haemophilus influenzae</i>	<i>Pseudomonas aeruginosa</i>
Human Metapneumovirus (hMPV)	<i>Legionella pneumophila</i>	<i>Staphylococcus epidermis</i>
Parainfluenza virus 1-4	<i>Mycobacterium tuberculosis</i>	<i>Streptococcus salivarius</i>

From *in-silico* analysis, it was concluded that no cross-reactivity is to be expected for the GSD NovaPrime® TSP SARS-CoV-2 primer and probe sequences.

In addition, wet testing was performed using the following clinical nasopharyngeal swab samples:

- 10 samples positive for RSV A
- 10 samples positive for RSV B
- 10 samples positive for human coronavirus 229E
- 20 samples positive for influenza B virus

No cross-reactivity was observed for the above-mentioned pathogens.

8.6. Matrix Comparison

8.6.1. Gargling Solution

Attention: Possible dilution effects and increased result variability when using gargling solution must be considered when evaluating the results. Likewise, the time of sampling (ideally directly after waking up without prior ingestion of liquids or food) and the duration of gargling with 4 mL of a 0.9 % NaCl solution (at least 30 seconds) are decisive for the result of the examination.

Method Comparison

The results obtained with the GSD NovaPrime® TSP SARS-CoV-2 were compared with the results of a CE marked reference RT-PCR assay.

		Reference (RT-PCR)	
		positive	negative
GSD NovaPrime® TSP SARS-CoV-2	positive	75	0
	negative	0	92

100 % agreement was achieved for all 167 samples tested.

LoD of NATtrol

The determinations were performed on the AriaDx™ (Agilent Technologies) Real-Time PCR instrument.

Negative gargling matrix was spiked with a dilution series of chemically inactivated intact SARS-CoV-2 virus particles (NATtrol™ SARS-CoV-2). Each dilution was extracted in triplicates using the GSD NovaPrime® IVD RNA Extraction AE1 kit. The provisional LoD was established by the lowest concentration that gives 100 % positive results when tested in triplicates. The final LoD was determined using a dilution series around the provisional LoD and was tested in 20 replicates in one run. The LoD is the lowest concentration demonstrating a positive signal in ≥ 95 % of 20 replicates tested.

Analyte	LoD**
NATtrol™ SARS-CoV-2	80 copies/mL

** Concentration values refer to the sample before RNA extraction

8.6.2. Saliva Solution

Attention: Sampling should be done more than 30 minutes after last ingestion of a beverage, food, chewing gum, cigarette/e-cigarette, tooth brushing or mouth rinse. After collection the samples need to be heat inactivated for 30 minutes at 56 °C, and subsequently 1 mL of PBS is added to the sample.

The high viscosity of the saliva matrix and variability in sample collection could result in high variance in viral load.

Method Comparison

59 positive and 59 negative saliva samples, as defined by a CE marked reference RT-PCR assay, were tested with the GSD NovaPrime® TSP SARS-CoV-2. The results obtained with the GSD NovaPrime® TSP SARS-CoV-2 were compared with the results of the reference RT-PCR assay.

		Reference (RT-PCR)	
		positive	negative
GSD NovaPrime® TSP SARS-CoV-2	positive	59	0
	negative	0	59

100 % agreement was achieved for all 118 samples tested.

Diagnostic Agreement

Matched pairs of nasopharyngeal swabs and saliva samples from 72 donors were tested using the GSD NovaPrime® RNA Extraction AE1 kit and the GSD NovaPrime® TSP SARS-CoV-2.

		Nasopharyngeal Swabs	
		positive	negative
Saliva	positive	44	4
	negative	5	19

Diagnostic Sensitivity: 89.8 % (95 % confidence interval: 81.6 % - 98.4 %)

Diagnostic Specificity: 82.6 % (95 % confidence interval: 67.7 % - 98.4 %)

LoD of NATtrol

The determinations were performed on the AriaDx™ (Agilent Technologies) Real-Time PCR instrument.

Negative saliva matrix was spiked with a dilution series of chemically inactivated intact SARS-CoV-2 virus particles (NATtrol™ SARS-CoV-2). Each dilution was extracted in triplicates using the GSD NovaPrime® IVD RNA Extraction AE1 kit. The provisional LoD was established by the lowest concentration that gives 100 % positive results when tested in triplicates. The final LoD was determined using a dilution series around the provisional LoD and was tested in 20 replicates in one run. The LoD is the lowest concentration demonstrating a positive signal in ≥ 95 % of 20 replicates tested.

Analyte	LoD**
NATtrol™ SARS-CoV-2	110 copies/mL

** Concentration values refer to the sample before RNA extraction

9. QUALITY CONTROL

In accordance with the manufacturer's ISO-certified Quality Management System, each lot of the GSD NovaPrime® TSP SARS-CoV-2 has been tested against predetermined specifications to ensure consistent product quality.

10. LIMITATIONS OF THE PROCEDURE

- Positive results indicate the presence of SARS-CoV-2, but do not rule out bacterial infection or co-infection with other pathogens not detected by the test.
- Very low levels of sample material may result in test runs with negative signals.
- Negative results do not exclude infection with SARS-CoV-2 and should not be used as the sole basis for treatment or other patient management decisions.
- Test results must be correlated with clinical observations, medical history, and/or epidemiologic information.

11. TRADEMARKS AND DISCLAIMERS

Registered names, trademarks, etc. used in this document are to be considered protected by law even if not specifically marked as such.

12. PRECAUTIONS AND WARNINGS

- The test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The use of the test kit with other analyzers than the ones mentioned in section 2. PRINCIPLE OF THE ASSAY has to be validated. Any change in design, composition, and test procedure as well as any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents for these reasons.
- Only for in-vitro diagnostic use.
- All patient samples should be regarded and handled as potentially infectious.
- Do not interchange reagents of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after expiry date stated on the label.
- Wear disposable powder-free gloves, a laboratory coat and eye protection when handling specimens.
- Always use DNase/RNase-free disposable reaction tubes and pipette tips with aerosol barriers.
- Avoid microbial and nuclease (DNase/RNase) contamination of the specimen and the components of the kit.
- In order to avoid contamination of working space with nucleic acids, reaction tubes/plates should not be opened after amplification.
- RT-PCR is highly sensitive to nucleic acid contamination. Therefore, positive / potentially positive material must be stored separately from all other components of the kit.
- Dedicate supplies and equipment to the separate working areas and do not move them from one area to another.
- This assay must not be used on the specimen directly.
- Prior to using this assay the nucleic acid has to be extracted with suitable extraction methods from the original specimen.
- Since ethanol is a strong Real-Time PCR inhibitor, it is necessary to completely eliminate it prior to the elution of the nucleic acid during extraction. If using spin columns with washing buffers containing ethanol, it is highly recommended to perform an additional centrifugation step of 10 min at approximately 17000 x g (~ 13000 rpm) before eluting the RNA. For this additional centrifugation step, use a new collection tube.
- Concentrations of > 4 % (v/v) ethanol resulted in inhibition of the GSD NovaPrime® TSP SARS-CoV-2 assay.
- The result of this RT-PCR kit may be influenced by potential mutations in the genome of the pathogen if they are located in the primer / probe binding region. Underestimation and/or failure to detect the pathogen may occur.
- RT-PCR inhibitors may also elicit underestimation, false negative results, or invalid runs. Therefore, only use nucleic acids extraction kits, which remove RT-PCR inhibitors, and which are dedicated for downstream RT-PCR processes.
- The Real-Time PCR is only designed for qualified personnel who are familiar with good laboratory practice and trained in RT-PCR.
- For further internal quality control each laboratory should additionally use known samples.

12.1. Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

For information about the packaging materials refer to PACKAGING MATERIALS.

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Gold Standard Diagnostics Frankfurt GmbH