17-OH Progesterone

Enzyme immunoassay for the quantitative determination of 17-OH Progesterone in human serum or plasma

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ENGLISH

1. INTRODUCTION

17-OH Progesterone (17-OH P) is an important intermediate of steroid biosynthesis. Elevated levels of 17-OH P are linked to dysfunction of steroid biosynthetic pathway that leads to an androgen excess. Therefore, measurements of 17-OH P are useful for the differential diagnosis of clinical conditions related to an hyperandrogenic phenotype and linked to adrenal disfunctions in both sexes.

Congenital Adrenal Hyperplasia (CAH) is a group of autosomal recessive disorders characterised by impaired cortisol synthesis¹⁻⁴. The predominant cause of the condition (95 % of cases) is mutations of the CYP21A2 gene that encodes the adrenal steroid 21-hydroxylase. 21-hydroxylase is responsible for conversion of 17-OHP to 11-deoxycortisol and progesterone to deoxycorticosterone. Assessment of 17-OH P levels provides supporting information for the differential diagnosis of CAH caused by mutations of the genes involved in cortisol synthesis from CYP21A2 gene mutations.

Measurement of 17-OH P is also considered useful in supporting diagnosis of other conditions characterised by hyperandrogenism such as polycystic ovarian syndrome (PCOS)⁵ and hirsutism in women, precocious puberty⁶ and precocious adrenarche. Determination of 17-OH P levels are used in such cases to exclude a diagnosis of CAH and non-classical CAH.

Patients affected by CAH shall be treated with glucocorticoids and mineral-corticoids. Treatment with glucocorticoids are of particular importance to avoid adrenal crisis and virilization due to excessive androgen synthesis^{2,4,7}. Since the goal of glucocorticoid treatments is not the complete suppression of 17-OH P synthesis, it is important to ensure appropriate monitoring of treatment to evaluate the correct dosage of the drugs. This is managed by the regular assessment of 17-OH P (and androstenedione) levels as traditional indicators of the adequacy of glucocorticoids treatment in CAH⁵.

2. INTENDED USE

Competitive immunoenzymatic colorimetric method (ELISA) for quantitative determination of 17-OH Progesterone in human serum or plasma (obtained by Lithium-heparine, Sodium-heparine or potassium EDTA).

3. PRINCIPLE OF THE ASSAY

Microtiter strip wells are precoated with anti-17-OH Progesterone antibodies (solid-phase). 17-OH Progesterone in the sample competes with added horseradish peroxidase labelled 17-OH Progesterone (enzyme-labelled antigen) for antibody binding. After incubation a bound/free separation is performed by solid-phase washing. The immune complex formed by enzyme-labelled antigen is visualized by adding Tetramethylbenzidine (TMB) substrate which gives a blue reaction product. The intensity of this product is inversely proportional to the amount of 17-OH Progesterone in the sample. Sulphuric acid is added to stop the reaction. This produces a yellow endpoint colour. Absorption at 450 nm is read using an ELISA microwell plate reader.

4. MATERIALS

4.1. Reagents supplied

- **Anti-17-OH Progesterone Coated Wells:** 12 break-apart 8-well snap-off strips coated with anti-17-OH Progesterone; in resealable aluminium foil.
- **Stop Solution:** 1 bottle containing 15 mL sulphuric acid, 0.15 mol/L (avoid any skin contact).
- **17-OH Progesterone-HRP Conjugate:** 1 bottle containing 22 mL of horseradish peroxidase labelled 17-OH Progesterone.
- **TMB Substrate Solution:** 1 bottle containing 15 mL 3,3',5,5'-tetramethylbenzidine (H₂O₂-TMB 0.26 g/L) (avoid any skin contact).
- **Wash solution 10x conc.:** 1 bottle containing 50 mL of a 10x concentrated solution of phosphate buffer 0.2 M, Proclin < 0.0015 %
- **17-OH Progesterone control A:** 1 bottle containing 1 mL of a lot-specific, ready to use control solution. The concentration is stated on the label.
- **17-OH Progesterone control B:** 1 bottle containing 1 mL of a lot-specific, ready to use control solution. The concentration is stated on the label.
- **17-OH Progesterone Standards:** 6 bottles, 1 mL each
 - Standard 0: 0 ng/mL
 - Standard 1: 0.2 ng/mL
 - Standard 2: 0.6 ng/mL
 - Standard 3: 2.0 ng/mL
 - Standard 4: 6.0 ng/mL
 - Standard 5: 16.0 ng/mL

4.2. Materials supplied

- 1 Cover foil
- 1 Instructions for use (IFU)

4.3. Materials and Equipment needed

- ELISA microwell plate reader, equipped for the measurement of absorbance at 450 nm, 620-630 nm
- Incubator 37 °C
- Manual or automatic equipment for rinsing wells
- Pipettes to deliver volumes between 10 and 1000 µL
- rotating mixer
- Distilled water
- Disposable tubes
- Timer

5. STABILITY AND STORAGE

Store all reagents at 2...8 °C in dark; the kit is stable up to the expiry date stated on label and on the Certificate of Analysis. Do not use the kit or its components after the expiry date.

Once opened, the kit is stable 6 months when stored at 2...8 °C.

6. REAGENT PREPARATION

It is very important to bring all reagents, samples and standards/controls to room temperature (22...28 °C) for at least 30 minutes before starting the test run! At the end of the assay, store immediately the reagents at 2...8 °C; avoid long exposure to room temperature.

6.1. Coated snap-off Strips

The ready to use break apart snap-off strips are coated with anti-17-OH Progesterone antibodies. Store at 2...8 °C. Open the bag only when it is at room temperature. Immediately after removal of strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C.

6.2. Conjugate

17-OH Progesterone-HRP Conjugate is a ready to use solution.

6.3. Standards

The standards are ready to use and have the following concentration of 17-OH Progesterone:

Standard 0:	0 ng/mL
Standard 1:	0.2 ng/mL
Standard 2:	0.6 ng/mL
Standard 3:	2.0 ng/mL
Standard 4:	6.0 ng/mL
Standard 5:	16.0 ng/mL

After first use the standards are still stable for another 6 months if stored at 2...8 °C, but no longer than the stated expiry date of the kit.

6.4. Controls

The bottles contain each 1 mL of a lot-specific control solution. The concentration is indicated on the label.

6.5. TMB Substrate Solution

The bottle contains 15 mL of a tetramethylbenzidine/hydrogen peroxide system. The reagent is ready to use and has to be stored at 2...8 °C in the dark. *The solution should be colourless or could have a slight blue tinge. If the substrate turns into blue, it may have become contaminated and should be thrown away.*

6.6. Stop Solution

The bottle contains 15 mL 0.15 M sulphuric acid solution. This ready to use solution has to be stored at 2...8°C.

6.7. Wash Solution

Dilute the content of the vial "10X Conc. Wash Solution" with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:10 dilution ratio. The diluted wash solution is stable for 30 days at 2...8 °C. In concentrated wash solution it is possible to observe the presence of crystals, in this case mix at room temperature until complete dissolution of crystals; for greater accuracy dilute the whole bottle of concentrated wash solution to 500 mL on taking care also transfer the crystals, then mix until crystals are completely dissolved.

7. SPECIMEN COLLECTION AND PREPARATION

Use human serum or plasma (obtained by Lithium-heparine, Sodium-heparine or potassium EDTA) samples with this assay. If the assay is performed within 4 days (96 hours) after sample collection, the samples should be kept at 2...8 °C; otherwise, they should be aliquoted and stored frozen. If samples are stored frozen, mix thawed samples well before testing. Avoid repeated freezing and thawing. Heat inactivation of samples is not recommended.

8. ASSAY PROCEDURE

8.1. Test Preparation

Please read the test instructions for use carefully **before** performing the assay. Result reliability depends on strict adherence to the instructions for use as described. Prior to commencing the assay, the distribution and identification plan for all specimens and standards/controls should be carefully established. Select the required number of microtiter strips or wells and insert them into the holder. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than one plate is used, it is recommended to repeat the dose response curve. Please allocate at least:

1 well	(e.g. A1)	for the substrate blank
2 wells	(e.g. B1+C1)	for standard 0
2 wells	(e.g. D1+E1)	for standard 1
2 wells	(e.g. F1+G1)	for standard 2
2 wells	(e.g. H1+A2)	for standard 3
2 wells	(e.g. B2+C2)	for standard 4
2 wells	(e.g. D2+E2)	for standard 5
2 wells	(e.g. F2+G2)	for control A
2 wells	(e.g. H2+A3)	for control B

It is necessary to determine standards/controls and patient samples in duplicate.

Perform all assay steps in the order given and without any appreciable delays between the steps.

A clean, disposable tip should be used for dispensing each standard/control and each patient sample.

1. Dispense 25 µL standards, controls and samples into their respective wells. Leave well A1 for substrate blank.
2. Add 200 µL 17-OH Progesterone-HRP Conjugate to each well. Leave well A1 for substrate blank.
3. Cover wells with the foil supplied in the kit.
4. **Incubate for 1 hour at 37 °C.**
5. When incubation has been completed, remove the foil, aspirate the content of the wells, and wash each well 6 times with 300 µL of diluted wash solution. Avoid overflows from the reaction wells. During each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.

Note: Washing is critical! Insufficient washing results in poor precision and falsely elevated absorbance values.

6. Dispense 100 µL TMB Substrate Solution into all wells.
7. **Incubate for exactly 15 min at room temperature (22...28 °C) in the dark.**
8. Dispense 100 µL Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate Solution. Shake the microplate gently.
Any blue colour developed during the incubation turns into yellow.
9. Measure the absorbance of the specimen at 450 nm against a reference wavelength of 620-630 nm or against blank within 5 minutes.

8.2. Measurement

Adjust the ELISA Microwell Plate Reader **to zero** using the **substrate blank in well A1**.

If - due to technical reasons - the ELISA reader cannot be adjusted to zero using the substrate blank in well A1, subtract the absorbance value of well A1 from all other absorbance values measured in order to obtain reliable results!

Measure the absorbance of all wells at **450 nm** and record the absorbance values for each standard and sample.

Where applicable calculate the **mean absorbance values** of all duplicates.

9. CALCULATION OF RESULTS

A variety of data reduction software packages are available, which may be employed to generate the mean calibration curve and to calculate the mean concentrations of unknown samples and controls. A 4-parameter logistic (4PL) curve fit, **including Standard 0 is required**.

Alternatively, a calibration curve may be prepared on semi-log graph paper by plotting mean absorbance on the Y-axis against concentration of analyte on the X-axis. Standard 0 should be included in the calibration curve. Read the mean absorbance value of each unknown sample off the curve.

10. QUALITY CONTROL

Good Laboratory Practice (GLP) requires the use of quality control specimens in each series of assays in order to check the performance of the assay. Controls should be treated as unknown samples, and the results analysed with appropriate statistical methods. The kit controls provided in the kit should be tested as unknowns and are intended to assist in assessing the validity of results obtained with each assay plate.

We recommend the users to maintain graphic records of the control values generated with each assay run, including the running means, SDs and %CVs. This information will facilitate the controls trending analysis relating to the performance of current and historical control lots. The trending will assist in the identification of assays which give control values significantly different from their average range.

When interpreting control data, users should note that this product was designed and developed as a manual product. The range stated on the Certificate of Analysis (CoA) should be appropriate for assays that are performed manually and with strict adherence to the Assay Procedure described above. It is recognised by Quality Control professionals, that as a result of differences in conditions and practices, there will always be variability in the mean values and precision of control measurements between different laboratories⁸.

11. MEASURING RANGE

The assay measuring range (AMR) is 0.26 – 16 ng/mL.

Any value that reads below 0.26 ng/mL should be reported as “< 0.26 ng/mL”. Any value that reads above 16 ng/mL should be reported as “> 16 ng/mL”.

12. METROLOGY AND TRACEABILITY

The 17-OH Progesterone ELISA has been standardised against internal reference standards (serum matrix) which have been value assigned to a commercially available LC-MS/MS method.

13. EXPECTED VALUES

The following ranges were determined using the 17-OH Progesterone ELISA and are provided for information only. The 95 % reference interval for apparently healthy adults were calculated by a non-parametric method following guidance from CLSI C28-A “Defining, Establishing and Verifying Reference Intervals in the Clinical Laboratory”.

		n°	Median (ng/mL)	Reference interval (ng/mL)
WOMEN	luteal phase	123	0.45	<LoQ – 2.19
	follicular phase	122	0.32	<LoQ – 2.10
	post-menopausal	122	<LoQ	<LoQ – 0.95
MEN		122	0.75	<LoQ – 1.97
CHILDREN (3 – 18 years)		124	<LoQ	<LoQ – 1.48

The above ranges should be considered as guidelines only; it is recommended that each laboratory establish its own expected range based upon its own patient population.

14. PERFORMANCE AND CHARACTERISTICS

14.1. Detection Capability

The limit of blank (LoB), limit of detection (LoD) and limit of quantitation (LoQ) were determined with guidance from CLSI EP17-A, “Protocols for Determination of Limits of Detection and Limits of Quantitation” using 6 blanks and 6 low level samples.

Sensitivity	Concentration
Limit of Blank (LoB)	0.02 ng/mL
Limit of Detection (LoD)	0.11 ng/mL
Limit of Quantitation (LoQ)	0.26 ng/mL

14.2. Trueness

Trueness has been demonstrated through method comparison of the 17-OH Progesterone ELISA to a commercially available LC-MS/MS using native donor samples – refer to section 14.5.

14.3. Precision

Precision of the 17-OH Progesterone ELISA was determined by performing a complex precision study.

Repeatability (within laboratory)

A total of 6 serum samples were assayed in 5 replicates, once a day for 5 days by 3 operators.

Data from one representative lot is shown below:

Sample	n	Mean Conc. (ng/mL)	Within run (Repeatability)	
			SD	CV%
1	75	0.90	0.10	11.4%
2	75	1.45	0.15	10.0%
3	75	2.10	0.20	9.3%
4	75	5.14	0.49	9.5%
5	75	8.94	0.74	8.3%
6	75	14.50	1.11	7.6%

Reproducibility

A total of 6 serum samples were assayed in 5 replicates, once a day for 5 days by 3 operators. Data from one representative lot is shown below:

Sample	n	Mean Conc. (ng/mL)	Within Laboratory (Reproducibility)	
			SD	CV%
1	75	0.90	0.12	13.0%
2	75	1.45	0.17	11.6%
3	75	2.10	0.23	11.1%
4	75	5.14	0.55	10.7%
5	75	8.94	0.82	9.2%
6	75	14.50	1.29	8.9%

14.4. Linearity

Linearity was evaluated based on CLSI EP-06, "Evaluation of the Linearity of Quantitative Measurement Procedures". For 17-OH Progesterone concentration by 17-OH Progesterone ELISA, the measurement procedure shows linearity for the interval from 0.2 to 17.22 ng/mL within the allowable deviation of linearity (ADL) of $\pm 15\%$.

14.5. Method comparison

The 17-OH Progesterone ELISA was compared against a commercially available LC-MS/MS method following CLSI EP-9A, "Method Comparison and Bias Estimation Using Patient Samples". A total of 120 samples, selected to represent a wide range of 17-OH Progesterone concentrations was assayed by each method. Passing-Bablok regression analysis was performed on the comparative data:

n	Slope [95% CI]	Intercept (ng/mL) [95% CI]	Correlation coefficient (r)
120	0.91 [0.87 to 0.97]	0.21 [0.16 to 0.28]	0.99

14.6. Cross-Reactivity

The analytical specificity was assessed with the following cross-reactants.

Cross-reactant	Concentration tested (unit)	Mean %Cross reactivity
11-deoxycortisol	80 ng/mL	2.15%
Progesterone	2000 ng/L	0.47%
Pregnenolone	2000 ng/mL	0.15%
Testosterone	2000 ng/mL	0.07%
17 β -Estradiol	2000 ng/mL	0.01%
Aldosterone	2000 ng/mL	0.00%
Estriol	2000 ng/mL	0.00%
Estrone-3-Sulphate	2000 ng/mL	0.00%
Spironolactone	2000 ng/mL	0.00%
Androstenedione	2000 ng/mL	0.05%
Androsterone	2000 ng/mL	0.04%
Corticosterone	2000 ng/mL	0.09%
Cortisol	2000 ng/mL	0.10%
Cortisone	2000 ng/mL	0.07%
DHEA	2000 ng/mL	0.07%
DHEA-S	20000 ng/mL	0.01%
DHT	2000 ng/mL	0.04%
Prednisolone	3000 ng/mL	0.02%
Prednisone	2000 ng/mL	0.01%

14.7. Interference

The following substances do not interfere with a bias of $\geq \pm 15\%$ in the 17-OH Progesterone ELISA when the concentrations are below the stated threshold presented in the following table.

Substance	Threshold Concentration
Bilirubin, conjugated	20 mg/dL
Bilirubin, unconjugated	20 mg/dL
Hemoglobin	200 mg/dL
Triglycerides	600 mg/dL

14.8. Serum-plasma study

The 17-OH Progesterone ELISA matrix comparison study was performed to evaluate the difference across tube types (serum separator tubes (SST), lithium heparin plasma, sodium heparin plasma and K2 EDTA plasma) versus the control samples (red top serum, without additive) following CLSI (EP9-A) guidelines. A total of 27 samples (23 native, 4 spiked) to cover the assay range were evaluated. Passing-Bablok regression analysis was performed on the comparative data:

Sample	Slope [95% CI]	Intercept (ng/mL) [95% CI]	Correlation coefficient (r)
SST	1.02 [0,91 to 1,13]	-0.03 [-0,21 to 0,11]	0.99
Lithium Heparin	1.00 [0,87 to 1,04]	0.00 [-0,07 to 0,08]	0.99
Sodium Heparin	1.00 [0,94 to 1,10]	-0.02 [-0,14 to 0,03]	0.99
EDTA	1.04 [0,90 to 1,12]	-0.03 [-0,19 to 0,09]	0.99

15. LIMITATIONS OF THE PROCEDURE

- As in the case of any diagnostic procedure, results must be interpreted in conjunction with the patient's clinical presentation and other information available to the physician.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays⁹. Patients routinely exposed to animals or to animal serum products can be prone to this interference and anomalous values may be observed.

16. PRECAUTIONS AND WARNINGS

- The test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The use of the testkits with analyzers and similar equipment has to be validated. Any change in design, composition and test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents for these reasons. The manufacturer is not liable for any results by visual analysis of the patient samples.
- Only for in-vitro diagnostic use by professional persons. Not for internal or external use in humans or animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- All components of human origin used for the production of these reagents have been tested for anti-HIV 1+2 antibodies, anti-HCV antibodies and HBsAg and have been found to be non-reactive. Nevertheless, all materials should still be regarded and handled as potentially infectious.
- Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not affected by BSE, but these materials should be handled as potentially infectious.
- Do not interchange reagents or strips of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after expiry date stated on the label.
- Use only clean pipette tips, dispensers, and lab ware.
- Do not interchange screw caps of reagent vials to avoid cross-contamination.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- After first opening and subsequent storage check conjugate and control vials for microbial contamination prior to further use.
- When using automatic equipment, the user has the responsibility to make sure that the kit has been appropriately tested. To improve the performance of the kit on ELISA automatic systems, it is recommended to increase the number of washes.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense conjugate without splashing accurately to the bottom of wells.
- Do not use heavily haemolysed or highly lipemic samples.
- Maximum precision is required for dispensation of the reagents.
- The TMB substrate contains an irritant, which may be harmful if inhaled, ingested, or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes. It may react with lead or copper plumbing to form potentially explosive metal azides. If you use a sink to remove the reagents, wash through with large amounts of water to prevent azide build-up.

- The stop solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of TMB substrate to direct sunlight, metal or oxidants. Do not freeze the solution.
- This method allows the determination of 17-OH Progesterone from 0.2 – 20 ng/mL.
- Treatment of the patient with cortisone, natural or synthetic steroids can impair 17-OH Progesterone determination.
- Some reagents contain small amounts of Sodium Azide or Proclin 300™ as preservative. Avoid the contact with skin or mucosa.
- Sodium Azide may be toxic if ingested or absorbed through the skin or eyes: moreover, it may react with lead or copper plumbing to form potentially explosive metal azides. If you use a sink to remove the reagents, allow scroll through large amounts of water to prevent azide build-up.
- The ELISA is only designed for qualified personnel following the standards of good laboratory practice (GLP).

WARNING: Sulphuric acid irritates eyes and skin. Keep out of the reach of children. Upon contact with the eyes, rinse thoroughly with water and consult a doctor!

16.1. Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

For information about the packaging materials refer to PACKAGING MATERIALS.

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
1. Honour JW. 17-Hydroxyprogesterone in children, adolescents and adults. *Ann Clin Biochem*. 2014 Jul;51(Pt 4):424-40.
2. Speiser PW, Arlt W, Auchus RJ, Baskin LS, Conway GS, Merke DP, Meyer-Bahlburg HFL, Miller WL, Murad MH, Oberfield SE, White PC. Congenital Adrenal Hyperplasia Due to Steroid 21-Hydroxylase Deficiency: An Endocrine Society Clinical Practice Guideline. *J Clin Endocrinol Metab*. 2018 Nov 1;103(11):4043-4088.
3. Bornstein SR, Allolio B, Arlt W, Barthel A, Don-Wauchope A, Hammer GD, Husebye ES, Merke DP, Murad MH, Stratakis CA, Torpy DJ. Diagnosis and Treatment of Primary Adrenal Insufficiency: An Endocrine Society Clinical Practice Guideline. *J Clin Endocrinol Metab*. 2016 Feb;101(2):364-89.
4. Sharma R, Seth A. Congenital adrenal hyperplasia: issues in diagnosis and treatment in children. *Indian J Pediatr*. 2014 Feb;81(2):178-85.
5. Legro RS, Arslanian SA, Ehrmann DA, Hoeger KM, Murad MH, Pasquali R, Welt CK; Endocrine Society. Diagnosis and treatment of polycystic ovary syndrome: an Endocrine Society clinical practice guideline. *J Clin Endocrinol Metab*. 2013 Dec;98(12):4565-92.
6. Grandone A, Marzuillo P, Luongo C, Toraldo R, Mariani M, Miraglia Del Giudice E, Perrone L. Basal levels of 17-hydroxyprogesterone can distinguish children with isolated precocious pubarche. *Pediatr Res*. 2018 Oct;84(4):533-536.
7. Rachoń D. Differential diagnosis of hyperandrogenism in women with polycystic ovary syndrome. *Exp Clin Endocrinol Diabetes*. 2012 Apr;120(4):205-9.
8. Basic QC Practices On-line Course; <http://www.Westgard.com>.
9. Boscatto, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27-33

PACKAGING MATERIALS/VERPACKUNGSMATERIALIEN/MATÉRIELS D'EMBALLAGE/MATERIALI D'IMBALLAGGIO/MATERIALES DE EMBALAJE/MATERIAIS DE EMBALAGEM

 PAP 21	 PAP 21	 PAP 22	MTP  ALU / LDPE 90
---	---	--	--

SYMBOLS KEY / SYMBOLSCHLÜSSEL / EXPLICATION DES SYMBOLES / LEGENDA / SÍMBOLOS / TABELA DE SÍMBOLOS

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17-OH Progesterone

Test Preparation

Prepare reagents and samples as described.
Establish the distribution and identification plan for all specimens and controls.
Select the required number of microtiter strips or wells and insert them into the holder.

Assay Procedure

	Substrate blank	Standard 0 - 5	Control A/B	Sample
Standard 0 - 5	-	25 µL	-	-
Control A/B	-	-	25 µL	-
Sample	-	-	-	25 µL
Conjugate	-	200 µL	200 µL	200 µL
Cover wells with foil supplied in the kit Incubate for 1 h at 37 °C Wash each well six times with 300 µL diluted wash solution. In case you automatic equipment, wash the wells at least 5 times				
TMB Substrate	100 µL	100 µL	100 µL	100 µL
Incubate for exactly 15 min at room temperature in the dark				
Stop Solution	100 µL	100 µL	100 µL	100 µL
Shake the microplate gently. Photometric measurement at 450 nm/620-630 nm				

Gold Standard Diagnostics Frankfurt GmbH

Waldstrasse 23 A6

63128 Dietzenbach, Germany

Tel.: +49 6074 23698-0

Fax: +49 6074 23698-900

E-Mail: info.frankfurt@eu.goldstandarddiagnostics.com

Website: clinical.goldstandarddiagnostics.com

DNOV004_17-OHP_IFU_rev01_fromLot_6009BN