
FREND Cardiac Triple

REF FRCT 025

- English
- German
- Italian
- Spanish
- French
- Portuguese

FREND Cardiac Triple

Quantitative Assay for Myoglobin, Troponin I, and
Creatine Kinase MB

REF FRCT 025

IVD For *in vitro* diagnostic use only

1. Intended use

The FREND Cardiac Triple test is a fluorescent nanoparticle immunoassay designed for *in vitro* quantitative measurement of myoglobin, troponin I, and creatine kinase MB in human plasma samples using EDTA as the anticoagulant. The FREND Cardiac Triple microfluidic flow cartridge is designed for use in the FREND™ System to aid in the diagnosis and assessment of the severity of myocardial infarction. FREND Cardiac Triple is for professional use only.

2. Summary and explanation of test

A leading cause of morbidity and mortality in the industrialized world is acute myocardial infarction (AMI).^[1] Emergency healthcare professionals regularly face the challenge of AMI diagnosis in patients with chest pain. Physical examination, patient history, electrocardiogram, and changes in serum biomarkers are key criteria for AMI diagnosis.^[2] In addition, rapid testing of serum cardiac markers to aid in guiding therapy has been shown to improve outcomes of high risk patients.^[3] When physical examination and other AMI indicators (like electrocardiogram, ECG) are inconclusive, levels of serum biomarkers such as the MB isoenzyme of creatine kinase (CK-MB), Myoglobin, and troponin I (TnI) can aid in helping a physician diagnose AMI.^[4]

CK-MB is an 82,000 dalton enzyme present in high concentrations in myocardial tissue.^[6] Elevated CK-MB serum concentration is commonly used to indicate AMI because of correlation to cardiac damage.^[6] 12 to 24 hours following AMI, CK-MB serum concentrations peak before returning to normal levels in approximately three days.^[7] Yet, because CK-MB elevation can result from events unrelated to cardiac muscle trauma, physicians rely on the concurrent measurement of additional cardiac biomarkers.^[8]

Cardiac troponin I (TnI) is a 24,000 dalton muscle regulation protein that becomes released into the blood stream upon AMI.^[9,10] Clinical data show elevation of cardiac TnI hours after AMI, but peaks at 8–24 hours before returning to normal levels days later.^[11] Cardiac TnI has high tissue specificity for AMI as well as other myocardial injuries such as unstable angina, cardiac contusions, cardiac transplant, coronary artery bypass graft surgery, congestive heart failure, or other physical trauma to the heart.^[12] While a variety of cardiac biomarkers are available for measurement, cardiac TnI has become the preferred biomarker for myocardial injury. Myocardial infarction has come to be defined as the rise and fall of cardiac biomarkers, such as TnI, in conjunction with ECG analysis and abnormal physical changes.^[13]

Myoglobin is a 17,800 dalton globular heme protein found in skeletal and cardiac muscle cells.^[14,15] Upon AMI, myoglobin is quickly released into the bloodstream, rapidly increasing in the first 2–4 hours, peaking at 8–10 hours, and returning to normal levels in 24 hours.^[16] Elevation of myoglobin, while quickly responsive to myocardial injuries, also has been observed in skeletal muscle damage, skeletal muscle disorders, renal failure, and strenuous exercise.^[17] Thus, the measurable multiple cardiac biomarkers can detect AMI in a broader time range despite varying tissue specificities.^[11]

3. Principle of the assay

The FRENDCardiac Triple test cartridge is a one-time-use rapid “sandwich” immunoassay utilizing fluorescent nanoparticles in microfluidic flow to capture and quantify myoglobin, TnI, and CK-MB in EDTA plasma specimens. A 35 µL drop of patient EDTA plasma sample is placed in the FRENDCardiac Triple cartridge, where the sample interacts with a proprietary mix of dry-loaded reagents. One of these reagents includes antibody-conjugated fluorescent nanoparticles, forming immune complexes with Myoglobin, TnI, and CK-MB in the patient sample. Capillary action moves the sample to the detection region, where capture antibodies grab the Myoglobin, TnI, or CK-MB nanoparticle complexes in detection zones. The fluorescence intensities of complexes are measured and Myoglobin, TnI, and CK-MB concentration is calculated by the FRENDCardiac™ System.

4. Material provided

Q'ty	Contents	Catalogue number
25	FREND Cardiac Triple Cartridge(s)	FRCT 025
30	Disposable pipette tip(s)	
1	FREND Cardiac Triple Code chip	
1	FREND Cardiac Triple Package Insert	

5. Materials required but not provided

- The FREND™ System
- Micro-pipette or any pipettor, capable of delivering 35 and 70 µL
- Personal protective equipment and biohazard waste disposal containers

6. Warning and Precautions

- For in vitro diagnostic use only.
- FREND Cardiac Triple cartridges are single-use.
- Use FREND Pipette and compatible tips, do not mouth pipette.
- Handle specimens in accordance with the OSHA Standard on Bloodborne Pathogens.^[18]
- Wash hands thoroughly after handling specimens.
- Wear protective gloves, clothing, and eyewear.
- Heterophilic Antibody Interference: Some individuals have antibodies to mouse protein (HAMA), goat protein (HAGA), or rabbit protein (HARA); interferences may occur in FREND Cardiac Triple.^[19,20]
- Do not use FREND Cardiac Triple cartridges beyond indicated expiration date.
- Dispose of used FREND Cardiac Triple cartridges and pipette tips as bio hazardous materials.
- Do not attempt to disassemble FREND Cardiac Triple cartridges: Micro-channels contain dry-loaded reagents that maybe bio-hazardous and/or toxic.
- Do not use FREND Cardiac Triple cartridge if pouch is damaged or seal is broken.

7. Storage and Stability

Store FREND Cardiac Triple cartridges at 2–8°C (35–46°F).

Do not remove the cartridge from sealed pouch until ready for use.

When stored properly, cartridges are stable until the marked expiration date.

Materials	Catalogue number
Refrigerator temperature storage (2~8 °C)	
FREND Cardiac Triple	FRCT 025
Room temperature storage (18~25 °C)	
Pipette tips	None

8. Specimen collection and handling

- Human EDTA plasma sample (including EDTA plasma collected in separator tubes) are suitable for use with FREND Cardiac Triple cartridges. Samples collected as serum in collection tubes or plasma in anticoagulant tubes (such as citrate or heparin) have not been validated for use in FREND Cardiac Triple cartridges.
- Follow instructions detailed in this package insert as well as the specimen collection tube manufacturer's instructions for specimen collection and preparation (including manufacturer's instructions for centrifugation time and speed).
- Collect blood by venipuncture technique into an EDTA plasma collection tube.
- After allowing the specimen to sufficiently mix with anticoagulant at room temperature, the sample tube can be centrifuged for 10 minutes at 3,000 rpm.
- The plasma can be immediately used or stored at -20 °C for 48 hours. Before performing a measurement on frozen sample using FREND Cardiac Triple cartridges, thaw the plasma specimen until equilibrated to room temperature. Avoid freeze-thaw cycles.
- For optimal results, avoid grossly hemolytic, lipidic, or turbid specimens.
- Specimens should be free of aggregated fibrin, red blood cells, or other particulate matter. When pipetting into the FREND Cardiac Triple cartridge sample inlet, ensure that bubbles in the sample are avoided. Bubbles may restrict flow and result in an incomplete or erroneous test result.

9. Procedure

Code chip installation

Each box of FREND Cardiac Triple cartridges contains a lot-specific FREND Cardiac Triple Code chip to correct for lot-to-lot deviation.

- 1) Press the 'Setup' button on the 'Main' screen.
- 2) Insert the FREND Cardiac Triple Code chip into the code chip slot, following the direction arrow, at the rear of the FREND™ System.
- 3) Press the 'Code chip' button on the 'Setup' screen.
- 4) The information from the FREND Cardiac Triple Code chip is automatically installed and saved.
- 5) When the installation is completed, press the 'OK' button to return to the 'Setup' screen,
- 6) Press the 'Item' button on the 'Setup'screen.
- 7) Confirm the FREND Cardiac Triple cartridge lot number and code chip installation date. Press the 'Home' button to return to the 'Main' screen.

Quality control

• FREND™ System QC cartridges

FREND™ QC Cartridge contains multiple controls to check optic part of the system. By testing QC Cartridge, part of analytical components of the system of (1) laser power, (2) alignment, and (3) mechanical integrity are confirmed.

The QC Cartridge should be run each day patient testing is planned. Refer to the Quality Control Procedures section in the User Manual of FREND™ System. In brief, perform QC Cartridge testing for the following conditions.

- 1) Upon initial setup of the system,
- 2) Each day of patient testing,
- 3) When the system has been transported or moved,
- 4) Whenever there is uncertainty about the performance of the system,
- 5) Whenever required by your laboratory's quality control requirements.

• Internal procedural controls

FREND Cardiac Triple cartridge contains built-in control features. Fluorescence signal in the Reference Zone of each cartridge shows: (1) that enough volume is added, (2) that proper flow is obtained, and (3) that the antibody is reactive. If this reference zone signal is missing or lower than threshold, the FREND™ System consider it as an incorrect or failed test, not

producing a test result but an error message. In addition, with each cartridge run, the system monitors, in part, for (1) flow of samples, (2) speed of sample flow, (3) shelf-life of cartridge components, (4) function of internal barcode scanner, and (5) function of scanner's mechanical components.

- **External quality control testing**

Commercially available controls from a variety of manufacturers are available that contain Cardiac Triple as a measured analyte. It is recommended that a minimum of at least two (2) levels of controls be run at least once per month or once for each new lot, whichever comes earlier. However, Controls should be run with a minimum frequency, depending on number of tests run in the laboratory. Each laboratory should establish its own criteria based on the following parameters:

- 1) Each new lot,
- 2) Each new shipment (even if from the same lot previously received),
- 3) Each new operator (an individual who has not run the tests for at least two weeks),
- 4) Monthly, as a continued check on storage conditions,
- 5) Whenever problems (storage, operator, or other) are identified,
- 6) Or other times as required by your laboratory's standard QC procedures.

Individual laboratory policy will dictate exactly which control materials and lot numbers should be run, the frequency with which controls are to be tested, criteria for acceptance of the results and required corrective action to be taken if results do not meet laboratory criteria. If any external quality control sample values are out of the acceptable range, it will be necessary to investigate the problem before reporting patient results to assure there is not an instrument or software malfunction. Do not assay patient samples on the FREND™ System using FREND Cardiac Triple if quality control results do not give expected values. Refer to your laboratory policies on how to determine acceptability of external control materials results. Each laboratory operates under a different set of regulations. Every laboratory must follow the standardized procedures acceptable to the regulatory agencies to which the laboratory is responsible.

Specimen processing

- **Preparation**

Remove from the refrigerator sufficient cartridges of FREND Cardiac Triple to test the number of patient samples and required external quality control materials. Allow the cartridges to come to room temperature for 15-30

minutes prior to the start of the testing sequence.

If using refrigerated patient samples, remove those from the refrigerator and allow them to come to room temperature prior to testing. If frozen samples will be utilized, be sure these are removed from the freezer, thawed naturally and then mixed gently but thoroughly prior to testing. Testing should not begin on these previously frozen samples until they have reached room temperature.

There are no other reagents or sample preparations necessary.

- **Assay procedure**

- 1) Prepare the FRENDCardiac Triple cartridge and specimen.
- 2) Record the Sample ID on the cartridge in the designated area.
- 3) Drop the sample (35 μ L) into sample inlet on the cartridges using a suitable micro-pipette equipped with a fresh pipette tip.
- 4) Press the 'Test' button on the 'Main' screen of the FRENDCardiac™ System.
- 5) The system moves to the Patient ID screen automatically.
- 6) Type the Patient ID and press the 'Enter' button to begin the test.
- 7) Insert the cartridge into the cartridge slot using the cartridge arrows as a guide.

Caution: Please check the direction of the cartridge before insertion and assure the insertion is complete.

- 8) When the reaction in the cartridge is complete, the FRENDCardiac™ System will automatically begin the reading process.
- 9) When the measurements are completed, the cartridge will automatically be expelled and the results displayed.

Caution: Do not remove power from the FRENDCardiac™ System while a cartridge is in the reading chamber. This may cause a system error.

- 10) If the FRENDCardiac™ System is connected to the optional printer, press the 'Print' button and the results will be output on the printer paper.
- 11) For more detailed instructions, please refer to the FRENDCardiac™ System User Manual

10. Procedural notes

If a Cardiac marker concentration(s) in a specimen is (are) found to be greater than the measurement range of 500 ng/mL (Myoglobin), 20 ng/mL (TnI), and 80 ng/mL (CK-MB), a dilution is (are) required.

The specimen should be diluted with low concentration sample that has been previously measured on the FREN D Cardiac Triple.

The recommended dilution ratio for the high and low concentration specimens is 1:10 (e.g. 10 μ L of high concentration specimen and 90 μ L of low concentration specimen).

It is desirable to dilute the sample so that the diluted specimen read between 50–200 ng/mL (Myoglobin), 1–10.00 ng/mL (TnI), and 5–50 ng/mL (CK-MB). Once a Cardiac marker result is obtained within the measurement range, the original concentration of the unknown can be easily manually calculated.

If a Cardiac marker concentration(s) in a specimen is(are) found to contain <5 ng/mL (Myoglobin), <0.05 ng/mL (TnI), or <1 ng/mL (CK-MB), then re-assay according to the 'Assay Procedure'.

※ Original unknown concentration of Cardiac marker = Concentration of diluted sample $\times$ 10.

11. Calculation of results

The FREN D™ System performs all samples and reagent handling operations automatically within the cartridge once the sample has been manually loaded to the sample inlet in the cartridge and the cartridge placed into the FREN D™ System. The rate of fluorescence produced by the reaction is read at various intervals during the analysis process, blank readings are subtracted after which the net rate is automatically converted to Cardiac Marker (Myoglobin, TnI and CK-MB) concentrations in ng/mL based upon information stored on the Cardiac Triple Code Chip. This result is then output on the screen and to the optional printer. It is also stored in memory on the FREN D™ System.

Screen displayed for various concentration scenarios

Displayed result	Description
	Myoglobin concentration Less than 5.00 ng/mL TnI concentration Less than 0.05 ng/mL CK-MB concentration Less than 1.00 ng/mL
	Myoglobin, TnI and CK-MB Concentration The result within the analytical measurement range
	Myoglobin concentration Higher than 500.00 ng/mL TnI concentration Higher than 20.00 ng/mL CK-MB concentration Higher than 80.00 ng/mL

12. Limitations of the procedure

- 1) The results from FRENDCARD Cardiac Triple cartridge should be used as a diagnostic aid to the clinician, in conjunction with the patient's medical history, other laboratory tests, and clinical examination.
- 2) Specimens can be measured with the reportable range of 5.00 ng/mL ~500.00 ng/mL myoglobin, 0.05 ng/mL ~ 20.00 ng/mL TnI, and 1.00 ng/mL ~80.00 ng/mL CK-MB. Samples below the reportable range will yield "<5.00 ng/mL myoglobin", "<0.05 ng/mL TnI", and/or "<1.00 ng/mL CK-MB" and samples above the reportable range will yield ">500.00 ng/mL Myoglobin", ">20.00 ng/mL TnI", and/or ">80.00 ng/mL CK-MB".

- 3) Specimens from patients with heterophilic antibodies, such as anti-mouse (HAMA), anti-goat (HAGA), or anti-rabbit (HARA) antibodies, maybe show falsely elevated or depressed values or may result in an incomplete test.^{11, 12} Patients routinely exposed to animals or animal serum products can be prone to these types of heterophilic interferences.
- 4) FREND Cardiac Triple is a sandwich immunoassay and can be affected by environmental conditions. Each laboratory must establish its own reference range based on its unique conditions and procedures.
- 5) CK-MB values may be elevated in patients with skeletal muscle injury as well as renal failure.
- 6) FREND Cardiac Triple cartridges have not been validated in point-of-care settings.

13. Expected value

EDTA plasma samples from a total of 276 normal, apparently healthy individuals (138 males and 138 females, Koreans) were assayed on 3 lots of the FREND Cardiac Triple assay using a single FREND™ System. The 99th percentiles are shown below:

Cardiac marker	Male (ng/mL)	Female (ng/mL)
Myoglobin	98.93	93.81
Troponin I	<0.05	<0.05
CK-MB	4.63	5.47

Results from individual laboratories may vary. Each laboratory should establish its own reference range based on the geographical, patient or environmental factor.

14. Performance characteristics

Precision

1) Within-run (Repeatability)

A precision study for FREND Cardiac Triple was performed as outlined in CLSI EP5-A2. Seven level standard samples were assayed in duplicates, twice per day over a period of 20 days.

1. Myoglobin

Sample ID	Conc. (ng/mL)	Within		Between-run		Between-day		Within-laboratory	
		SD	CV(%)	SD	CV(%)	SD	CV(%)	SD	CV(%)
A	0.0	0.258	46.2	0.065	11.6	0.081	14.5	0.278	49.8
B	26.63	2.265	8.6	1.066	4.1	0.872	3.3	2.651	10.1
C	53.13	3.151	5.9	1.563	2.9	1.163	2.2	3.704	7.0
D	86.25	4.530	5.2	0.437	0.5	1.860	2.1	4.916	5.7
E	152.50	6.824	4.5	0.141	0.1	3.861	2.5	7.842	5.2
F	285.00	13.698	4.8	2.162	0.8	4.868	1.7	14.697	5.1
G	550.00	28.232	5.1	6.045	1.1	0.712	0.1	28.880	5.2

2. Troponin I

Sample ID	Conc. (ng/mL)	Within		Between-run		Between-day		Within-laboratory	
		SD	CV(%)	SD	CV(%)	SD	CV(%)	SD	CV(%)
A	0.00	0.006	65.0	0.001	11.6	0.002	18.0	0.007	68.5
B	0.30	0.029	9.4	0.008	2.6	0.005	1.5	0.030	9.8
C	1.50	0.086	5.7	0.006	0.4	0.007	0.5	0.086	5.7
D	3.00	0.167	5.5	0.026	0.9	0.045	1.5	0.175	5.8
E	6.00	0.273	4.5	0.241	4.0	0.150	2.5	0.394	6.6
F	12.00	0.607	5.0	0.147	1.2	0.171	1.4	0.647	5.3
G	24.00	1.115	4.6	0.138	0.6	0.282	1.2	1.158	4.8

3. CK-MB

Sample ID	Conc. (ng/mL)	Within		Between-run		Between-day		Within-laboratory	
		SD	CV(%)	SD	CV(%)	SD	CV(%)	SD	CV(%)
A	0.00	0.104	49.5	0.024	11.3	0.034	16.4	0.112	53.3
B	1.10	0.103	9.5	0.044	4.1	0.023	2.1	0.115	10.6
C	5.50	0.296	5.4	0.121	2.2	0.123	2.3	0.343	6.3
D	11.00	0.743	6.8	0.383	3.5	0.098	0.9	0.842	7.7
E	22.00	1.265	5.7	0.147	0.7	0.513	2.3	1.373	6.2
F	44.00	2.746	6.2	0.565	1.3	0.964	2.2	2.964	6.7
G	88.00	4.562	5.1	1.737	2.0	1.913	2.2	5.243	5.9

2) Reproducibility (Lot)

The results met with acceptance criteria in reproducibility test of FREND Cardiac Triple.

1. Myoglobin

Sample ID	Conc. (ng/mL)	Mean concentration (ng/mL)			Mean	SD	CV(%)
		Lot 152003	Lot 152004	Lot 152005			
A	0.0	0.31	0.32	0.53	0.38	0.122	31.7
B	26.63	26.76	28.32	27.54	27.54	0.780	2.8
C	53.13	53.43	54.39	53.01	53.61	0.706	1.3
D	86.25	85.88	91.91	82.74	86.84	4.660	5.4
E	152.50	152.85	161.88	143.95	152.89	8.965	5.9
F	285.00	283.26	274.94	301.53	286.58	13.602	4.7
G	550.00	557.91	535.79	585.62	559.77	24.966	4.5

2. Troponin I

Sample ID	Conc. (ng/mL)	Mean concentration (ng/mL)			Mean	SD	CV(%)
		Lot 152003	Lot 152004	Lot 152005			
A	0.00	0.01	0.01	0.01	0.01	0.001	23.1
B	0.30	0.30	0.31	0.31	0.31	0.008	2.6
C	1.50	1.46	1.48	1.45	1.46	0.018	1.2
D	3.00	2.92	2.98	2.92	2.93	0.035	1.2
E	6.00	5.88	5.85	5.97	5.92	0.063	1.1
F	12.00	12.23	12.44	12.01	12.18	0.215	1.8
G	24.00	22.96	23.65	24.58	23.94	0.813	3.4


3. CK-MB

Sample ID	Conc. (ng/mL)	Mean concentration (ng/mL)			Mean	SD	CV(%)
		Lot 152003	Lot 152004	Lot 152005			
A	0.00	0.26	0.17	0.21	0.21	0.045	21.4
B	1.10	1.06	1.11	1.12	1.10	0.032	2.9
C	5.50	5.54	5.48	5.26	5.38	0.143	2.7
D	11.00	10.61	10.87	11.35	11.05	0.375	3.4
E	22.00	22.51	21.72	23.20	22.66	0.739	3.3
F	44.00	41.66	44.87	45.16	44.21	1.942	4.4
G	88.00	85.23	89.59	82.89	85.15	3.401	4.0


Dilution linearity

The dilution linearity study as outlined in CLSI EP6-A was performed by diluting a high concentration Cardiac Triple specimen with a low concentration Cardiac Triple specimen. Results are summarized in the graphs below, but may vary in individual laboratories.


1) Myoglobin

2) Troponin I

3) CK-MB

The analytical measurement range of FREN D Cardiac Triple, Myoglobin is from 5.00 ng/mL to 500.00 ng/mL, Troponin I is from 0.05 ng/mL to 20.00 ng/mL, CK-MB is from 1.00 ng/mL to 80.00 ng/mL.

Analytical sensitivity

The limit of blank (LoB) and limit of detection (LoD) was determined using guidelines found in CLSI document EP17-A. LoB was determined using 12 replicate measurements of five blank patient samples (0.00 ng/mL). LoD was determined using 12 replicates measurements of five low level patient samples.


Cardiac marker	LoB (ng/mL)	LoD (ng/mL)	LoQ (ng/mL)
Myoglobin	1.72	3.14	4.30
Troponin I	0.028	0.040	0.041
CK-MB	0.01	0.26	0.41

Comparative analysis


FREND Cardiac Triple was compared to the Predicate device Cardiac Triple using guidelines outlined in CLSI guideline EP9-A2. Clinical samples (EDTA plasma, n=94) were measured in duplicate on both systems. There are high correlation value ($0.95 \leq r$) between FREND Cardiac Triple and Predicate device.

Cardiac marker	Slope	Correlation coefficient (r)
Myoglobin	1.0634	0.9771
Troponin I	0.9044	0.9761
CK-MB	0.9284	0.9739


1) Myoglobin

2) Troponin I

3) CK-MB

Specificity

The following substances were evaluated for potential cross-reactivity with FRENDCardiac Triple at the concentration indicated below. Testing was done according to the instruction recommended in CLSI guideline EP7-A2. No significant cross-reactivity was found.

Substance	Concentration (ng/mL)
Actin	1,000
CK-BB	1,000
CK-MM	1,000
Cardiac Troponin C	2,000
Cardiac Troponin T	2,000
Myosin	2,000
Skeletal Troponin I	1,000
Skeletal Troponin T	1,000
Tropomyosin	2,000

Interference

The interference study was performed as recommended in the CLSI guideline EP7-A2 using four concentrations of Myoglobin, Troponin I, CK-MB. No interference by the substances below was found. A variety of common medications were tested for interference with the FREND Cardiac Triple. The testing showed there was no significant interference from the tested drugs that would affect the interpretation of Myoglobin, Troponin I, CK-MB result as assayed on the FREND Cardiac Triple.

No.	Substances	Concentration
1	Hemoglobin	1,000 mg/dL
2	Bilirubin	20 mg/dL
3	Triglycerides	1,000 mg/dL
4	Albumin	5,000 mg/dL
5	Acetaminophen Bioextra	1,324 µmol/L
6	Allopurinol	294 µmol/L
7	Atenolol	37.6 µmol/L
8	Caffeine	308 µmol/L
9	Captopril	23 µmol/L
10	Chloramphenicol	155 µmol/L
11	Diclofenac sodium salt	169 µmol/L
12	Digoxin	7.8 µmol/L

13	Diltiazem hydrochloride	15 µmol/L
14	Enalapril maleate salt	0.86 µmol/L
15	Erythromycin	81.6 µmol/L
16	Hydrochlorothiazide	20.2 µmol/L
17	Indomethacin	100 µmol/L
18	Nifedipine	1,156 µmol/L
19	Probenecid	2,100 µmol/L
20	Propamolol hydrochloride	7.71 µmol/L
21	Quinidine	37 µmol/L
22	Sulfamethoxazole	1.12 µmol/L
23	Theophylline	222 µmol/L
24	Trimethoprim	138 µmol/L
25	Verapamil hydrochloride	4.4 µmol/L

15. References

- 1) American Heart Association: Cardiovascular disease statistics.
- 2) Ryan, T.J., Antman, E.M., Brooks, N.H., et al. ACC/AHA guidelines for the management of patients with acute myocardial infarction: A report of the American College of Cardiology/American Heart Association Task Force on Practice Guidelines (Committee on Management of Acute Myocardial Infarction). Journal of the American College of Cardiology. 1999, 34, 890-911.
- 3) Apple, F.S., Christenson, R.H., Valdes, R., et al. Simultaneous rapid measurement of whole blood myoglobin, creatine kinase MB, and cardiac troponin I by the triage cardiac panel for detection of myocardial infarction. Clinical Chemistry. 1999, 45:2, 199-205.
- 4) Morrow, D.A. Cannon, C.P., Jesse, R.L., et al. National Academy of Clinical Biochemistry Laboratory Medicine Practice Guidelines: Clinical Characteristics and Utilization of Biochemical Markers in Acute Coronary Syndromes. Clinical Chemistry. 2007, 53:4, 552-574.
- 5) Neumeier, D. Tissue specific and subcellular distribution of creatine kinase isoenzymes. Creatine Kinase Isoenzymes. New York: Springer-Verlag, 1981, 85-131.
- 6) Tsung, S.H. Creatine kinase isoenzyme patterns in human tissue obtained at surgery. Clinical Chemistry. 1976, 22/2, 173-175.
- 7) Lott, J.A. Serum enzyme determinations in the diagnosis of acute myocardial infarction: and update. Human Pathology. 1984, 15:8, 706-716.
- 8) Morrow, D.A., Cannon, C.P. Chest pain and biomarkers: has troponin supplanted CK-MB, ACC Current Journal Review. 2001, 10, 18-21.


- 9) Mair, J., Wanger, I., Pushendorf, B., et al. Cardiac troponin I to diagnose myocardial injury. *The Lancet*. 1993, 341:8848, 838-839.
- 10) Bertinchant, J.P., Larue, C., Pernel, I., et al. Release kinetics of serum cardiac troponin I in ischemic myocardial injury. *Clinical Biochemistry*. 1996, 29:6, 587-594.
- 11) Mair, J., Morandell, D., Genser, N., et al. Equivalent early sensitivities of myoglobin, creatine kinase MB mass, creatine kinase isoform ratios, and cardiac troponin I and T for acute myocardial infarction. *Clinical Chemistry*. 1995, 41:9, 1266-1272.
- 12) Mair, J., Genser, N., Morandell, D., et al. Cardiac troponin I in the diagnosis of myocardial injury and infarction. *Clinica Chimica Acta*. 1996, 245, 19-38.
- 13) National Academy of Clinical Biochemistry. Laboratory Medicine Practice Guidelines: Clinical Characteristics and Utilization of Biochemical Markers in Acute Coronary Syndromes. *Clinical Chemistry*. 2007, 53:4, 552-574.
- 14) Vaidya, H.C. Myoglobin. *Lab Medicine*. 1992, 23, 306-310.
- 15) Bhayana, V., Henderson, A.R. Biochemical markers of Myocardial Damage. *Clinical Biochemistry*. 1995, 28:1, 1-29.
- 16) Rozenman, Y., Gotsman, M. The earliest diagnosis of acute myocardial infarction. *Annual Review of Medicine*. 1994, 45, 31-44.
- 17) Bhayana, V., Cohoe, S., Pellar, T.G., et al. Combination testing for myocardial infarction using myoglobin, CK-2 (mass), and troponin I. *Clinical Biochemistry*. 1994, 27:5, 395-406.
- 18) 29 CFR Part 1910.1030 Bloodborne Pathogens. US Department of Labor, Occupational Safety and Health Administration.
- 19) Schroff, R.W. et al. Human Anti-murine immunoglobulin Responses in Patients Receiving Monoclonal Antibody Therapy. *Cancer Research*. 1985, 45, 879-885.
- 20) Boscatto, L.M. et al. Heterophilic antibodies: A Problem for All Immunoassay. *Clinical Chemistry*. 1988, 34:1, 27-33.

Glossary of Symbols

	Caution, warning, Consult accompanying documents
	Catalogue number/Reference number
 www.nanoentek.com/eifu.php	Consult Instructions for Use An electronic instructions for use (eIFU) indicator (website address) may accompany the symbol when used to indicate an instruction to consult an eIFU.
	Lot number/Batch number
	Use by YYYY-MM-DD or YYYY-MM
	Manufacturer
	CE marking
	<i>In vitro</i> diagnostic medical device
	Temperature limitation
	Contains sufficient for <n> tests
	Do not reuse
	Do not use if package is damaged
	For prescription use only CAUTION: Federal (U.S.) law restricts this device to sale by or on order of a physician.
	US Corporation
	Patient ID
	Result
	Sample Drop
	Authorized representative in the European Community

UK Representative	Authorized representative in United Kingdom
CH REP	Authorized representative in Switzerland
BRH	Authorized representative in Brazil

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15. Referenzen

- 1) American Heart Association: Cardiovascular disease statistics.
- 2) Ryan, T.J., Antman, E.M., Brooks, N.H., et al. ACC/AHA guidelines for the management of patients with acute myocardial infarction: A report of the American College of Cardiology/American Heart Association Task Force on Practice Guidelines (Committee on Management of Acute Myocardial Infarction). Journal of the American College of Cardiology. 1999, 34, 890-911.
- 3) Apple, F.S., Christenson, R.H., Valdes, R., et al. Simultaneous rapid measurement of whole blood myoglobin, creatine kinase MB, and cardiac troponin I by the triage cardiac panel for detection of myocardial infarction. Clinical Chemistry. 1999, 45:2, 199-205.
- 4) Morrow, D.A. Cannon, C.P., Jesse, R.L., et al. National Academy of Clinical Biochemistry Laboratory Medicine Practice Guidelines: Clinical Characteristics and Utilization of Biochemical Markers in Acute Coronary Syndromes. Clinical Chemistry. 2007, 53:4, 552-574.
- 5) Neumeier, D. Tissue specific and subcellular distribution of creatine kinase isoenzymes. Creatine Kinase Isoenzymes. New York: Springer-Verlag, 1981, 85-131.

- 6) Tsung, S.H. Creatine kinase isoenzyme patterns in human tissue obtained at surgery. *Clinical Chemistry*. 1976, 22/2, 173-175.
- 7) Lott, J.A. Serum enzyme determinations in the diagnosis of acute myocardial infarction: and update. *Human Pathology*. 1984, 15:8, 706-716.
- 8) Morrow, D.A., Cannon, C.P. Chest pain and biomarkers: has troponin supplanted CK-MB, ACC Current Journal Review. 2001, 10, 18-21.
- 9) Mair, J., Wanger, I., Pushendorf, B., et al. Cardiac troponin I to diagnose myocardial injury. *The Lancet*. 1993, 341:8848, 838-839.
- 10) Bertinchant, J.P., Larue, C., Pernel, I., et al. Release kinetics of serum cardiac troponin I in ischemic myocardial injury. *Clinical Biochemistry*. 1996, 29:6, 587-594.
- 11) Mair, J., Morandell, D., Genser, N., et al. Equivalent early sensitivities of myoglobin, creatine kinase MB mass, creatine kinase isoform ratios, and cardiac troponin I and T for acute myocardial infarction. *Clinical Chemistry*. 1995, 41:9, 1266-1272.
- 12) Mair, J., Genser, N., Morandell, D., et al. Cardiac troponin I in the diagnosis of myocardial injury and infarction. *Clinica Chimica Acta*. 1996, 245, 19-38.
- 13) National Academy of Clinical Biochemistry. Laboratory Medicine Practice Guidelines: Clinical Characteristics and Utilization of Biochemical Markers in Acute Coronary Syndromes. *Clinical Chemistry*. 2007, 53:4, 552-574.
- 14) Vaidya, H.C. Myoglobin. *Lab Medicine*. 1992, 23, 306-310.
- 15) Bhayana, V., Henderson, A.R. Biochemical markers of Myocardial Damage. *Clinical Biochemistry*. 1995, 28:1, 1-29.
- 16) Rozenman, Y., Gotsman, M. The earliest diagnosis of acute myocardial infarction. *Annual Review of Medicine*. 1994, 45, 31-44.
- 17) Bhayana, V., Cohoe, S., Pellar, T.G., et al. Combination testing for myocardial infarction using myoglobin, CK-2 (mass), and troponin I. *Clinical Biochemistry*. 1994, 27:5, 395-406.
- 18) 29 CFR Part 1910.1030 Bloodborne Pathogens. US Department of Labor, Occupational Safety and Health Administration.
- 19) Schroff, R.W. et al. Human Anti-murine immunoglobulin Responses in Patients Receiving Monoclonal Antibody Therapy. *Cancer Research*. 1985, 45, 879-885.
- 20) Boscatto, L.M. et al. Heterophilic antibodies: A Problem for All Immunoassay. *Clinical Chemistry*. 1988, 34:1, 27-33.

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15. Riferimenti

- 1) American Heart Association: Cardiovascular disease statistics.
- 2) Ryan, T.J., Antman, E.M., Brooks, N.H., et al. ACC/AHA guidelines for the management of patients with acute myocardial infarction: A report of the American College of Cardiology/American Heart Association Task Force on Practice Guidelines (Committee on Management of Acute Myocardial Infarction). Journal of the American College of Cardiology. 1999, 34, 890-911.
- 3) Apple, F.S., Christenson, R.H., Valdes, R., et al. Simultaneous rapid measurement of whole blood myoglobin, creatine kinase MB, and cardiac troponin I by the triage cardiac panel for detection of myocardial infarction. Clinical Chemistry. 1999, 45:2, 199-205.
- 4) Morrow, D.A. Cannon, C.P., Jesse, R.L., et al. National Academy of Clinical Biochemistry Laboratory Medicine Practice Guidelines: Clinical Characteristics and Utilization of Biochemical Markers in Acute Coronary Syndromes. Clinical Chemistry. 2007, 53:4, 552-574.

- 5) Neumeier, D. Tissue specific and subcellular distribution of creatine kinase isoenzymes. *Creatine Kinase Isoenzymes*. New York: Springer-Verlag, 1981, 85-131.
- 6) Tsung, S.H. Creatine kinase isoenzyme patterns in human tissue obtained at surgery. *Clinical Chemistry*. 1976, 22/2, 173-175.
- 7) Lott, J.A. Serum enzyme determinations in the diagnosis of acute myocardial infarction: and update. *Human Pathology*. 1984, 15:8, 706-716.
- 8) Morrow, D.A., Cannon, C.P. Chest pain and biomarkers: has troponin supplanted CK-MB, *ACC Current Journal Review*. 2001, 10, 18-21.
- 9) Mair, J., Wanger, I., Pushendorf, B., et al. Cardiac troponin I to diagnose myocardial injury. *The Lancet*. 1993, 341:8848, 838-839.
- 10) Bertinchant, J.P., Larue, C., Pernel, I., et al. Release kinetics of serum cardiac troponin I in ischemic myocardial injury. *Clinical Biochemistry*. 1996, 29:6, 587-594.
- 11) Mair, J., Morandell, D., Genser, N., et al. Equivalent early sensitivities of myoglobin, creatine kinase MB mass, creatine kinase isoform ratios, and cardiac troponin I and T for acute myocardial infarction. *Clinical Chemistry*. 1995, 41:9, 1266-1272.
- 12) Mair, J., Genser, N., Morandell, D., et al. Cardiac troponin I in the diagnosis of myocardial injury and infarction. *Clinica Chimica Acta*. 1996, 245, 19-38.
- 13) National Academy of Clinical Biochemistry. Laboratory Medicine Practice Guidelines: Clinical Characteristics and Utilization of Biochemical Markers in Acute Coronary Syndromes. *Clinical Chemistry*. 2007, 53:4, 552-574.
- 14) Vaidya, H.C. Myoglobin. *Lab Medicine*. 1992, 23, 306-310.
- 15) Bhayana, V., Henderson, A.R. Biochemical markers of Myocardial Damage. *Clinical Biochemistry*. 1995, 28:1, 1-29.
- 16) Rozenman, Y., Gotsman, M. The earliest diagnosis of acute myocardial infarction. *Annual Review of Medicine*. 1994, 45, 31-44.
- 17) Bhayana, V., Cohoe, S., Pellar, T.G., et al. Combination testing for myocardial infarction using myoglobin, CK-2 (mass), and troponin I. *Clinical Biochemistry*. 1994, 27:5, 395-406.
- 18) 29 CFR Part 1910.1030 Bloodborne Pathogens. US Department of Labor, Occupational Safety and Health Administration.
- 19) Schroff, R.W. et al. Human Anti-murine immunoglobulin Responses in Patients Receiving Monoclonal Antibody Therapy. *Cancer Research*. 1985, 45, 879-885.
- 20) Boscato, L.M. et al. Heterophilic antibodies: A Problem for All Immunoassay. *Clinical Chemistry*. 1988, 34:1, 27-33.

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NanoEntek, Inc.
MT Promedt Consulting GmbH