





## References

1. Lars-Olof Hanson et al. Current Opinion In Infectious Diseases 1997; 10:196-201.
2. Chetana Vaishnavi. Immunology and Infectious Diseases 1996; 6:139-144.
3. Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACC Press, 1995.
4. Yoshitsugu Hokama et al. Journal of Clinical Lab. Status 1987; 1:15-27.
5. Kari Pulki et al. Scand J Clin Lab Invest 1986; 46:606-607.
6. Werner Muller et al. Journal of Immunological Methods 1985; 80:77-90.
7. Shogo Otsuji et al. Clin Chem 1982; 28/10:2121-2124. Clinical Chemistry

AR. No.: AR I84

ASLCRP1-2023-04-001

ISO 13485 CERTIFIED COMPANY 

 Accurex Biomedical Pvt. Ltd.

Plant - G-54/55, MIDC Tarapur, Boisar, Palghar - 401 506. Maharashtra, INDIA

E-mail: help@accurex.net

Website: www.accurex.net

Call us toll free on : 74004 53000

Code: <b>ACRP</b>	Category: <b>Open</b>	Group: <b>ACCUREX 2</b>
Name: <b>ACRP</b>	Type: <b>Fixed time</b>	

  

<b>Calibration</b>	<b>Reagent blank</b>	<b>Correlation</b>	<b>Settings</b>	<b>Edit</b>			
<b>General</b>		<b>Reagent blank</b>	<b>Validity</b> <input type="text" value="360"/> hours	<b>Blank correction</b> <input type="checkbox"/>			
Number of standards: <input type="text" value="1"/>							
Validity (hours): <input type="text" value="360"/>							
<b>Standards</b>	<b>Single</b>	<b>Calibrator</b>	<b>Dilution rate</b>	<b>Lot</b>	<b>Primary concnt</b>	<b>Final concnt</b>	<b>Reps.</b>
1	CRPCAL		1	20200930	60.0000	60.0000	2

<b>ACRP</b>	<b>TIME</b>	<b>WASHING</b>	<b>AUTODILUTION</b>
Category: <b>Open</b>	Incubation 1: <input type="text" value="145"/> sec	Needle: <input type="text" value="3"/>	Enable autodilution: <input type="checkbox"/>
Double click category to edit	Incubation 2: <input type="text" value="29"/> sec	Cuvette: <input type="text" value="3"/>	Rate: <input type="text" value="4"/>
Version date: <b>4/3/2023 3:26:01 PM</b>	Reading: <input type="text" value="116"/> sec	Use cleaning solution: <input type="checkbox"/>	Max OD: <input type="text" value="0"/> Abs
Type: <b>Fixed time</b>	<b>DILUTION</b>	Needle post-wash: <input type="text" value="0"/>	Max Conc: <input type="text" value="150"/> mg / L
Name: <b>ACRP</b>	Serum	<b>LIMITS</b>	Min Conc: <input type="text" value="0"/> mg / L
Version: <b>1</b>	<input checked="" type="checkbox"/> 1:1 <input type="checkbox"/> 1:2 <input type="checkbox"/> 1:4	Blank OD min: <input type="text" value="0.4"/> Abs	<b>REFLEX</b>
Group: <b>ACCUREX 2</b>	<input type="checkbox"/> 1:10 <input type="checkbox"/> 1:100 <input type="checkbox"/> 1:40	Blank OD max: <input type="text" value="0.9"/> Abs	Enable reflex: <input type="checkbox"/>
<b>REAGENTS</b>	Custom rate: <input type="text" value="40"/>	Reaction slope: <b>Positive</b>	Method: <input type="text"/>
Reagent 1 ID: <b>CRP1</b>	QC dilution: <input type="text" value="1"/>	Calibration slope: <b>Positive</b>	Limit type: <b>Deactivated</b>
Reag 1 bottle: <b>Large</b>	<b>Urine</b>	OD Range min: <input type="text" value="0"/> Abs	Conc. limit: <input type="text" value="0"/> mg / L
Reagent 2 ID: <b>CRP2</b>	<input type="checkbox"/> 1:1 <input type="checkbox"/> 1:2 <input type="checkbox"/> 1:4	OD Range max: <input type="text" value="2.5"/> Abs	<b>OUTLIERS</b>
Reag 2 bottle: <b>Small</b>	<input type="checkbox"/> 1:10 <input type="checkbox"/> 1:100 <input type="checkbox"/> 1:40	Blank OD Δ max: <input type="text" value="0"/> Abs	Max. curve deviation: <input type="text" value="5"/> %
<b>PIPETTING</b>	Custom rate: <input type="text" value="40"/>	Min conc: <input type="text" value="0"/> mg / L	<b>PATHOLOGICAL RANGES</b>
Type: <b>RL + R2 / S</b>	QC dilution: <input type="text" value="1"/>	Max conc: <input type="text" value="150"/> mg / L	Mini... <input type="text" value="0"/> mg / L
Sample vol: <input type="text" value="2"/> μL	Diluent vol: <input type="text" value="5"/> μL	<b>UNITS</b>	Sample type: <input type="text"/>
Reagent 1 vol: <input type="text" value="240"/> μL	<b>WAVELENGTHS</b>	Units 1: <input type="text" value="mg / L"/>	Maximum... <input type="text" value="6"/> mg / L
Reagent 2 vol: <input type="text" value="60"/> μL	Wavelength 1: <input type="text" value="600"/> nm	Units 2: <input type="text"/>	
		Units 3: <input type="text"/>	
		Decimal digits: <input type="text" value="2"/>	
		Target molecule: <input type="text"/>	
		Display as inhibition: <input type="checkbox"/>	
			<b>Add Remove Edit</b>