


Introduction


1. AutoPure ALT (GPT) is reagent kit for direct quantitative determination of alanine transaminase (ALT) in human serum and plasma on automated clinical chemistry analyzers.
2. AutoPure ALT (GPT) is a ready-to-use, two liquid reagent system.
3. With AutoPure ALT (GPT), the assay is linear upto 600 IU/l (10.02 µkat/l).

Principle

α-Ketoglutarate reacts with L-alanine in the presence of ALT (GPT) to form pyruvate and L-glutamate.



The increase in pyruvate is determined in an indicator reaction catalyzed by lactate dehydrogenase.



The conversion of NADH to NAD⁺ is directly proportional to the ALT (GPT) activity in the specimen and is measured kinetically.

Reagent Storage, Stability & Handling

AutoPure ALT (GPT) is a ready-to-use, two liquid reagent system.

Shelf life

Stable till the expiry date indicated on the label, when stored at 2°–8°C.

On-Board Reagent Stability

R1: 30 days at 2°–8°C after opening.

R2: 30 days at 2°–8°C after opening.

Protect the reagent from light and contamination.

Do not freeze the reagent.

Components & Concentration of Working Solution


Component	Concentration
R1	
• Tris buffer	100 mmol/l
• L-alanine	≥ 500 mmol/l
• Lactate dehydrogenase	≥ 3000 IU/l
• Stabilizers, excipients & surface active agents.	
R2	
• Tris buffer	20 mmol/l
• NADH	≥ 0.18 mmol/l
• α-Ketoglutarate	≥ 10 mmol/l
• Stabilizers, excipients & surface active agents.	

Specimen Collection & Preservation

Collect sample using standard sampling tube. Although serum is preferred, heparinized plasma or EDTA plasma can be used. Specimen with any visible haemolysis is not acceptable. ALT (GPT) activity in serum/plasma is stable for 7 days at 2°–8°C and 30 days at -10°C. Centrifuge samples containing precipitate before performing the assay.

Procedure

AutoPure ALT (GPT) can be used on various automated analyzers. The procedure described below is for Sphera auto-analyzer.



Calculations

Fully automated system automatically calculates the ALT (GPT) activity of each sample.

Results in µkat/l = Results in IU/l × 0.0167

Application Sheet

Refer to the application sheet for details. For additional system applications, contact our local Accurex representative.

Calibration

For calibration, it is recommended to use C.f.a.s* or any other suitable calibrator material

Calibration frequency

Re-calibration is recommended

- Whenever the reagent lot is changed.
- As per the requirement of quality control procedures.

*C.f.a.s is brand name of Roche Diagnostics, Germany

Quality Control

Each batch of AutoPure ALT (GPT) is assayed with atleast six quality control sera – Precinorm**, Precipath**, Biorad I***, Biorad II***, Accutrol normal and Accutrol Abnormal prior to release*.

To ensure the adequate quality control, it is recommended that the laboratory should use the normal and abnormal commercial reference control serum. It should be realized that the use of quality control material checks both the reagent and instrument functions together.

If the control values fall outside the specified limits, each of the below criteria should be cross-checked and corrected:

- Proper instrument function – wavelength setting, light source and temperature control.
- Cleanliness of probes and cuvettes.
- Bacterial contamination of wash water used by the instrument.
- Expiry date of the reagent kit.

Brand name of Roche Diagnostics, Germany; *Brand name of Biorad, USA; *subject to availability of control serum.

Expected Values

Serum/ Plasma

Men : Upto 41 IU/l (0.68 µkat/l)

Women : Upto 31 IU/l (0.52 µkat/l)

Note:

Expected range varies from population to population. It is therefore recommended that each laboratory should establish its own range. For diagnostic purposes, the ALT (GPT) results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

Performance Characteristics

Linearity

With AutoPure ALT (GPT), the assay is linear upto 600 IU/l (10.02 µkat/l). Determine samples with higher concentrations via the rerun function. On instruments without rerun function, manually dilute samples with higher concentrations using 0.9% NaCl or distilled/deionized water (e.g. 1 + 3). Multiply the result by the appropriate dilution factor (e.g. 4).

Interference

There is no significant interference in samples containing upto 60 mg/dl of bilirubin. Avoid haemolysis as it interferes with the assay.

Precision

Reproducibility was determined using quality control sera as shown below:

n=15

Quality Control Material	Within run			Between run		
	Mean IU/l	SD IU/l	%CV	Mean IU/l	SD IU/l	%CV
Low Control Serum	32.7	0.55	1.7	32.5	0.49	1.5
High Control Serum	115	1.27	1.1	115	1.05	0.9

Co-Relation Studies

A comparison of the ALT (GPT) determination using AutoPure ALT (GPT) and AutoZyme SGPT reagent gave the following co-relation (IU/l):

Linear Regression

y = 0.994x + 0.383

r = 0.9942

Sy.x = 0.86

No. of samples measured : 65

The sample concentrations measured were between 11 and 40 IU/l.

References

1. Tietz N.W., ed. *Clinical Guide to Laboratory Tests*, 3rd ed. Philadelphia, Pa: W.B. Saunders, 1995: 20-21.

Application sheet



- Bergmeyer H.U., Horder M, Rej R. Approved recommendation. (1985) on IFCC methods for the measurement of catalytical concentration of enzymes, Part 3. IFCC method for L- alanine aminotransferase. *J.Clin. Chem. Clin.Biochem* 1986; 24: 481-489.
- Fischbach F, Zawta B. Age – dependent Reference Limits of Serval Enzymes in Plasma at Different Measuring Temperatures. *Clin Lab.* 1992; 38 : 555 – 561.
- Penttila , I.M. et al., Scand, *J. Clin.Lab. Invest* 35, 275 (1975).
- Hafkensheid, J. C.M. et al. *J. Clin. Chem. Clin. Biochem.* 17, 219 (1979).

Designations																		
Code: ALT	Name: GPT/ALT	Type: Kinetic	Group: ACCUREX															
1 - Pipetting Reagent 1 ID: ALT R1 Reagent 1 bottle: Large Reagent 2 ID: ALT R2 Reagent 2 bottle: Small Sample vol: 1st <input type="text"/> μ l 2nd <input type="text"/> μ l Reagent 1 vol: 200 μ l Reagent 2 vol: 50 μ l Diluent vol: <input type="text"/> μ l	3 - Wavelengths Wavelength 1: 340 nm Wavelength 2: 0 nm 4 - Washing Needle: 1 Cuvette: 1 5 - Incompatibility 1 <input type="text"/> <input type="text"/> 2 <input type="text"/> <input type="text"/> 3 <input type="text"/> <input type="text"/> 4 <input type="text"/> <input type="text"/>	7 - Autodilution Rate: * <input type="text"/> <input type="text"/> Max OD: ** <input type="text"/> <input type="text"/> Abs 8 - Dilutions Serum <input checked="" type="checkbox"/> 1:1 <input type="checkbox"/> 1:2 <input type="checkbox"/> 1:4 <input type="checkbox"/> 1:10 <input type="checkbox"/> 1:40 <input type="checkbox"/> 1:100 Urine <input checked="" type="checkbox"/> 1:1 <input type="checkbox"/> 1:2 <input type="checkbox"/> 1:4 <input type="checkbox"/> 1:10 <input type="checkbox"/> 1:40 <input type="checkbox"/> 1:100	9 - Pathological ranges <table border="1"> <thead> <tr> <th>Minimum</th> <th>Sample type</th> <th>Maximum</th> </tr> </thead> <tbody> <tr> <td></td> <td>Male</td> <td>41</td> </tr> <tr> <td></td> <td>Female</td> <td>31</td> </tr> <tr> <td></td> <td></td> <td></td> </tr> <tr> <td></td> <td></td> <td></td> </tr> </tbody> </table> Add Remove Edit 10 - Results units Units 1: <input type="text"/> IU/L Units 2: <input type="text"/> Conversion : <input type="text"/> 1 Decimal digits: <input type="text"/> 1	Minimum	Sample type	Maximum		Male	41		Female	31						
Minimum	Sample type	Maximum																
	Male	41																
	Female	31																
2 - Time Incubation 1: 145 Sec Incubation 2: 58 Sec Reading: 87 Sec	6 - Links Blank OD min: <input type="text"/> 0.7 Abs Blank OD max: <input type="text"/> 2 Abs Reaction slope: Negative OD Range min: <input type="text"/> 0 Abs OD Range max: ** 0.17 Abs	Linearity: <input type="text"/> 25 % Subst Depletion: * <input type="text"/> 0.2 Abs Blank OD: <input type="text"/> 0 Abs Min conc: <input type="text"/> 0 IU/L Max conc: <input type="text"/> 600 IU/L																

* User defined ** User defined based on calibration

Designations			
Code: ALT	Name: GPT/ALT	Type: Kinetic	Group: ACCUREX
OK			
Calibration			
General		Number of standards: 3	Factor: <input type="text"/>
		Validity (hours): 240	Interpolation type: Average
Standards		Current Calibration	
<input type="button" value="Single"/>	Calibrator	Dilution rate	Lot
	Primary concent	Final concent	OD value
	Date		
1	CFAS	<input type="text"/> 1	1 <input type="text"/>
2	CFAS	<input type="text"/> 1	2 <input type="text"/>
3	CFAS	<input type="text"/> 1	3 <input type="text"/>
Reagent blank <input type="checkbox"/>	Blank Correction <input type="checkbox"/>	Validity 9999 hours	OD <input type="text"/> OD 0 <input type="text"/> Date <input type="text"/>
Correlation			
Y = <input type="text"/> 1.000 X + <input type="text"/> 0.000			
<input type="button" value="Print calibration"/>		<input type="button" value="Curve chart"/>	

AR. No.: I42A

ASLPI-2014-11-001

ISO 13485, ISO 9001 CERTIFIED COMPANY

Accurex Biomedical Pvt. Ltd.

Head Office - Mumbai. Tel.: 91 (022) 67446744; Fax: 91 (022) 67446755
 E-mail: accurex@vsnl.com; Website: www.accurex.org
 Plant : G-54, MIDC Tarapur, Boisar, Thane - 401 506. INDIA.
 Call us toll free on : 1800 209 8456