





Biocard® HCV RAPID TEST KIT

Biocard® HCV Rapid is qualitative in-vitro diagnostic test for the detection of antibodies to Hepatitis C Virus



Approved by National Institute Biologicals (NIB)



Sensitivity : 100% Specificity : 100%



Presence of control line to validate the assay



Sample types: Serum, Plasma (10µL)/ Whole Blood (20µL)



Pre coated with HCV Core NS3, NS4 and NS5



Results with in 20 minutes




Storage: 4°C to 30°C



Biocard®


HCV RAPID TEST KIT

Test Procedure



Results Interpretation

For Biocard® HCV



*refer kit insert for exact procedure

Performance

Status	Positive	Negative
Positive	110	0
Negative	0	650
Total	110	650

Diagnostic Sensitivity: 100 %

Diagnostic Specificity: 100 %

Ordering Information

Product Code	Product Name	Pack Size	Storage	Shelf Life
30201416	Biocard® HCV RAPID TEST KIT	30 Test	4-30°C	24 Months

*Please note: Product specifications are subject to change without prior notice owing to product modifications, improvements / up-gradation.