

COVIDsure Pro

MULTIPLEX RT PCR KIT

- Simultaneous detection of 3 specific COVID-19 genes ORF1ab, N & E besides internal control
- Low virus load detection: < 5 copies/ μ l
- Extracted RNA samples from entire respiratory tract can be used
- Compatible with multiple RT-PCR Instruments

Approved by



icmr
INDIAN COUNCIL OF
MEDICAL RESEARCH
Serving the nation since 1911

NIV
NATIONAL INSTITUTE
OF VIROLOGY

COVIDsure PRO MULTIPLEX RT PCR KIT

COVIDsure Pro is a 4 Gene Kit. It is developed for the rapid detection of viral infection by simultaneous targeting of three major genes ORF 1ab, envelope (E), and nucleocapsid (N). The Kit also includes an internal control to ensure successful extraction and amplification of Viral RNA.

SPECIFICATIONS



Highly Specific

- No detectable cross reactivity with most human Coronavirus strains, respiratory pathogens, etc
- Specificity - 100%

Utmost Sensitivity

- Wider range of Covid-19 strains covered to detect most mutations (mutation variants)
- Sensitivity - 100%



Time-Saving Approach

- Processing time is approximately 1 hour

Certified Quality

- CE IVD marked reliable quality
- Approved & Validated by ICMR



Ordering Information

Product Code	Product Description	Pack Size	Storage	Shelf Life
30201297	COVIDsure Pro Multiplex RT PCR Kit	50 Tests	-20°C	12 months
30201298	COVIDsure Pro Multiplex RT PCR Kit	100 Tests	-20°C	12 months



*Please note: Product specifications are subject to change without prior notice owing to product modifications, improvements / up-gradation.

corporate@trivitron.com | www.trivitron.com

Customer Care Helpline: +91 98400 80008



TRIVITRON
HEALTHCARE
speaking your language