


# Thrombospan LS aPTT - EA

Ellagic acid activated reagent for Partial Thromboplastin time



## Reliable & Accurate Results

- Rabbit Brain Phospholipid Cephaloplastin reagent ensures ease in end point detection by manual method and minimal interference by photo optical instruments
- Sensitive and responsive to Lupus Anticoagulant to prevent it (if present) in interfering with factor assays
- Additional 3.2 % Sodium Citrate which is recommended by WHO (World Health Organization) as anticoagulant for blood collection along with the kit

## High Sensitivity

- Ellagic Acid as a Contact Activator, ensures efficient activation of factor XII, quick initiation of intrinsic pathway
- Sensitive to heparin

## Convenient and Time Saving

- Liquid-stable, ready to use two reagents
- High Reagent Shelflife upto 15 Months at +2 °C to +8 °C

## Flexibility

- Suitable for manual procedure as well as on semi automated and fully automatated coagulation analysers
- Pack Size of 6 mL (3 X 3 mL), suitable for low and medium workload laboratory

## Product Ordering Details:

Product Code	Pack-Size	Shelf Life
98LS200-06	3 X 3 mL	15 Months

For commercial and product related query, please contact :

**arkray healthcare Pvt. Ltd.**

701/702, Opulence, 6th Road, TPS III, Santacruz (East), Mumbai - 400 055.

Telephone: +91-22-6155 9355, Fascimile: +91-22-6155 9356, URL <http://www.arkray.co.in>

Specification may change without prior notice

ISO 9001:2008 / ISO 13485:2003 / GMP