

# *Acculyte Touch*

## *Acculyte-3P & 5P*



*a life saving partner*

Aspen Diagnostics Pvt. Ltd.  
**(RAPID DIAGNOSTIC GROUP OF COMPANIES)**  
ISO 9001 : 2015 CERTIFIED COMPANY

# Acculyte Touch



Acculyte-3P



Acculyte-5P

## Features:

- Intuitive graphical user interface, touch screen and graphically guided instructions
- Fast and error-free data entry using barcode scanner
- Sample size: 65µl, and fast test results: 30 seconds
- Big volume reagent pack is optional
- Low analyzer purchase price and low cost per test
- Automatic or on-demand calibration
- Sleep mode to save reagent
- High efficiency cleaning procedure
- On-line operation and self-diagnostic tutorials
- Above 1000 results storage



## Sample Type

Plasma, Serum, Whole Blood, Urine

## Models

K<sup>+</sup>/Na<sup>+</sup>/Cl<sup>-</sup>    K<sup>+</sup>/Na<sup>+</sup>/Cl<sup>-</sup>/Ca<sup>2+</sup>/pH

## Test Range and Accuracy

| Analyte          | Measure Range (mmol/l) | Resolution (mmol/l) | Precision C.V. |
|------------------|------------------------|---------------------|----------------|
| Na <sup>+</sup>  | 20.0~200.0 mmol/L      | 0.1                 | <1.0%          |
| K <sup>+</sup>   | 0.50~10.00 mmol/L      | 0.01                | <1.0%          |
| Cl <sup>-</sup>  | 20.0~200.0 mmol/L      | 0.1                 | <1.0%          |
| Ca <sup>2+</sup> | 0.3~5.0 mmol/L         | 0.1                 | <1.5%          |
| pH               | 6.00~9.00              | 0.01                | <1.0%          |

## Working Environment

|                     |                       |
|---------------------|-----------------------|
| Relative Humidity   | ≤ 85%                 |
| Ambient Temperature | 50C~400C              |
| Sun Radiation       | Avoid direct sunshine |
| Power Supply        | 100-240V ~ 50/60Hz    |

\* The product specifications are subject to change without prior notice owing to product modification, information/upgradation



Imp & Mktg.  
Aspen Diagnostics Pvt. Ltd.  
**(RAPID DIAGNOSTIC GROUP OF COMPANIES)**  
**ISO 9001 : 2015 CERTIFIED COMPANY**  
Mfd.: M.C. Hi-Tech Co., Ltd., P.R. Ch.

Local Contact:

|  |
|--|
|  |
|--|