




# ***Stat Test for Malaria...***



*Early*  
**Detection**

*Correct*  
**Diagnosis**

*Right*  
**Treatment**



# ***FalciVax***


---

***....fulfill the needs***

# FalciVax

## Rapid Test for Malaria Antigen Pv/Pf


RDT for *P. vivax* malaria & *P. falciparum* malaria with direct speciation.




- True speciation of the “Big Two” i.e. *P. falciparum* & *P. vivax* - through specific detection of Pf HRP - 2 & specific detection of Pv pLDH.
- Detection of mixed infection of *P. falciparum* & *P. vivax*.
- Sensitivity and specificity in excess of 98% for *P. falciparum* & *P. vivax*.
- No interference of Heterophile antibodies in test results.
- Well suited for field conditions, on-site diagnosis, pathological labs and reference labs, alike - uses whole blood from finger prick or primary tube, as a sample.
- Easy to use- rapid, 20 minutes assay.



### INTERPRETATION OF RESULT



Test negative for malaria



Test positive for *P. f.*



Test positive for *P. v.*



Test positive for mixed infection



Test invalid

Malaria  
RDT's that

SEEK • DETECT • ENSURE

ISO 13485  
Certified Company



T TULIP DIAGNOSTICS (P) LTD.

For further information contact :

Gitanjali, Tulip Block, Dr. Antonio do Rego Bagh, Alto Santacruz, Bambolim Complex Post Office, Goa - 403 202, INDIA.  
Tel.: +91 832 2458546-50 Fax : +91 832 2458544 E-mail : sales@tulipgroup.com Website : www.tulipgroup.com