

# XPRESS GLUCO+

BLOOD GLUCOSE MONITORING SYSTEM



## Reliable Diabetes Monitoring

-  Results in 5 Seconds
-  Avg reading upto 90 days
-  0.5µl Sample Volume

# XPRESS GLUCO

BLOOD GLUCOSE MONITORING SYSTEM



TREND  
INDICATOR



GOD  
METHOD



AUTO  
DETECTION  
OF STRIP



EXCELLENT  
PRECISION &  
ACCURACY

**“From #1 GLUCOSE lab test REAGENT manufacturer Since 1984”**

SPECIFICATIONS	
Principle	Electrochemical Biosensor
Automatic Shutoff	Yes
Battery Life	1000 Tests
Hematocrit Range	15% - 65%
Measuring Units	mg/dL & mmol/L
Sample Volume	0.5 $\mu$ L
Sample Type	Whole Blood
Shelf Life of Strips	24 months (From manufacture date)
Available pack Size of Strips	50 & 100
Result Range	20-800 mg/dL
HCT Interference	No
Memory	1000 Tests
Operating Temperature	2°C-50°C
Storage Condition	Room Temperature

## HEAD OFFICE

212 Udyog Mandir No. 1, 7-C Bhagoji Keer Marg Mahim (w),  
Mumbai - 400016. Maharashtra, India

 +91 7400453000

 [help@accurex.net](mailto:help@accurex.net)

## BRANCH OFFICE

86-90 Paul Street, London,  
England, EC2A 4NE



  
Since 1984