

uRIQUICK Pro

Wireless | Compact | Advanced

Introducing
the next-generation
Urine Analyzer



**2.8" Colour
Touch Screen**



**Wireless
LIS**



**Auto
Calibration**



**Closed
System**



**Built In
Thermal Printer**

Step into the Future with Wireless Connectivity

uRIQUICK Pro Urine Analyzer

Testing Range

Up to
14 Parameters



RFID Based
Closed
System

Albumin-
to-creatinine
ACR ratio

Throughput
Up to
120 Tests/Hr



2.8" Colour
Touch Screen



WiFi
Enabled

HL7

Wireless LIS
As per Standard



Urine Analyzer Technical Characteristics Comparison

FEATURES	STRUMed URIQUICK Pro	EUROPEAN ANALYZERS	AMERICAN ANALYZERS	CHINESE ANALYZERS
Reflectance Photometry Principle				
Complete Visible Spectrum Detection (400nm - 700nm)				
Wireless LIS				
Colour Display				
Unlimited Memory Storage via File Dump				
Internal Barcode Scanner				
Local Spare Parts Availability				
Replacement Warranty				
Low Cost Maintenance				
Indigenous Manufacturer				

Uriquick Pro Urine Strips -
100 Strips / Container

PARAMETERS	Leukoocytes
	Nitrites
	Urobilinogen
	Protein
	pH
	Blood
	Specific Gravity
	Ketone
	Bilirubin
	Glucose
	Microalbumin
	Ascorbate
	Creatinine
	Calcium