


TRUSTgen HCV Real Time PCR Kit (Quantitative)



CODE#	PACK SIZE	SPECIMEN	SHELF LIFE
AP1010	25 Test	Serum/Plasma	12 Months

Hepatitis C Virus (HCV) is a blood-borne virus that can cause both acute and chronic liver infections, potentially leading to cirrhosis, liver failure, or liver cancer. TRUSTgen HCV Real Time PCR Quantitative Kit is developed with enhanced limit of detection and higher linearity, aiding in early diagnosis and monitoring treatment response.

TRUSTgen HCV Real Time PCR Kit (Quantitative)

- Includes detection of HCV – 1a, 1b, 2b, 3a, 4a, 5a and 6a genotypes
- Limit of Detection: 10 IU/mL (6 copies/mL)
- Linearity range: 40 IU/mL to 1E+08 IU/mL
- Compatible with all the PCR instruments & extraction systems

FEATURES



Ready to use
quantitative standards



Shorter time to
results



Sensitivity – 100%
Specificity – 100%

