

Smart healthcare, BUDDI™ checks your patient

BUDDI™



Bright field **U**tilization **D**irect **D**iagnostic Immunoassay analyzer

- **Clarifying decision making**
- **Accurate and fast result**
- **Smart connectivity**




Test Item

- **Infectious Disease**
Influenza A & B
Strep A
RSV
- **Cardiac disease**
H-FABP / Troponin I Duo
- **Asthma**
EDN (Eosinophil Derived Neurotoxin)
- **Gastrointestinal infection**
Adenovirus / Norovirus / Rotavirus

BUDDI™ is portable solution that is able to read results of various diseases with rapid test kit.



- » Provide objective criteria
- » User friendly touch screen interface
- » Store over 3,000 tests
- » Easy data management

Smart connectivity



User friendly interface

Offers Intuitive & easily Accessible GUI



Specification

Light source	Display	Interface	Hardware platform		Operating system	Battery
Green LED module	7" TFT LCD	Touch screen	Tablet PC		Windows 8.1	Up-to 4 hours (Built in battery)
Electrical rating	Channel	Connectivity	Test ID	Camera	Dimension	Weight
AC100 ~ 240V, 50~60Hz	Single	Wi-Fi, SMS, email	2D Barcode	CMOS	166 mm x 202 mm X 205 mm (WxLxH)	2.0 kg