



ISO 9001:2015
QUALITY MANAGEMENT SYSTEM

ISO 13485:2016
MEDICAL DEVICE QUALITY MANAGEMENT

DENGUE NS1 RAPID TEST KIT



Fast



Accurate



Precise




info@promea.in

www.promea.in



DENGUE NS1 RAPID TEST KIT

3 Steps Quick Procedure



Technical Specifications

- Lateral Flow Immuno Chromatography Assay.
- Detects NS1 Antigen.
- Quick qualitative results.
- Easy storage and handling.
- Operating Temperature 2°C to 30°C.
- Simple, convenient and easy to interpret.
- High Precision and Accuracy.
- Target Customers Primary Health Centre, Government Public Health Programs, Pharma Industries, Research and Other Affected Countries.
- Specificity 96.15%.
- Sensitivity 98.3%.
- Designed against Indian Population.
- Serum/Plasma Samples.
- Sample Size 60µL.



✉ bdm@promea.in | ☎ 91549 38441

🏭 Manufactured by:

Promea Therapeutics Pvt Ltd.

Plot No. M4, Medical Device Park, Sultanpur (V), Patancheru (M),
Sangareddy (Dist), Telangana - 502319. www.promea.in

promea®