




CoviDx™ mPlex-4R

SARS-CoV-2 RT-PCR Detection Kit

Multiplex Detection & Identification of 3 COVID-19 Target genes E, N, RdRP & 1 Human internal Control RNaseP



Features



Highly Efficient

PCR with high specificity & sensitivity




Fastest Process

Results within 80 minutes post-extraction



Multiplex Real Time PCR

Detection & Identification of 3 COVID Target Gene & 1 Control Gene in a single tube




Mighty Mastermix

PCR results Highest inhibitor tolerance, Thermostability & Reproducibility



Varied Specimens

Human Respiratory Sample Specimen



Simple workflow

User Friendly convenient workflow

Order information



Cat No

CM4R01



Product Name

CoviDx Mplex-4R SARS-Cov-2-RT-PCR Detection Kit



RxN

100

Manufactured by:

NeoDx Biotech Labs Pvt. Ltd. 24, NH 44, Pragati Plaza, 4th Floor, opp Biocon, Veer Sandra, Hebbagodi, Bengaluru, Karnataka 560100

Distributed By:

DSS Imagetech Pvt. Ltd. A-5, Mohan Co-op Industrial Estate, Mathura Rd, New Delhi, 110044

info@neodx.in
www.neodx.in