


iNSTA XPLOR P - PA001

SEMI AUTOMATIC PROTEIN ANALYZER



Precision Proteins, Delivered Simply




**OPTO
MECHANICAL
TECHNOLOGY**



AUTO CALIBRATION



**INTERNAL
INCUBATION
SYSTEM**



**PID CONTROLLED
HEATING SYSTEM**

STRUMed Solutions Private Limited
Plot No. 44, Ichhapore GIDC Industrial Estate, Bhatpore,
Surat, Gujarat - 394510 | Contact:- +91 6379454154
Email:- customercare@strumed.com | Web:- www.strumed.com



iNSTA XPLOR P - PA001

SEMI AUTOMATIC PROTEIN ANALYZER



FEATURES	SPECIFICATIONS
PRINCIPLE	TURBO-DENSITOMETRY AND NEPHELOMETRY WITH OPTO-MECHANICAL SYSTEM
TEST PARAMETERS	CRP, FERRITIN, HbA1c, MICROALBUMIN, RF
TEST SELECTION	*UPCOMING - D DIMER, CYSTATIN C, ASO, hsCRP, PCT, VITAMIN D, ANTI CCP, IgG, IgA, IgM, IgE
TESTING CHANNELS	RFID CARD BASED SELECTION WITH CLOSED SYSTEM
MEMORY	500 PATIENT RECORDS
LIS CONNECTIVITY	WI-FI BASED WIRELESS LIS SUPPORT
BARCODE SCANNER	EXTERNAL BARCODE SCANNER
PRINTER	INTERNAL PRINTER – AUTOMATIC RESULT PRINT AVAILABLE
POWER SUPPLY	INPUT: 150-240V, OUTPUT: 12V-3A
DISPLAY	2.8" TFT TOUCH SCREEN COLOR DISPLAY
DIMENSION	300 X 300 X 165 MM (L X B X H)
WEIGHT	3 KG

