


Universal Pipette Tips



Pipette tips Volume Filter/Non filter

10ul -39mm Pipette tips Filter/Non filter
 10ul -47mm XL Pipette tips Filter
 20ul -56mm Pipette tips Filter
 100ul -56mm Pipette tips Filter
 200ul -56mm Pipette tips Filter
 1000ul -78mm Pipette tips Filter
 1250ul -100mm Pipette tips Filter

STERILE FILTER TIPS, RACKED, IN FILP LID BOX

Filtered Tip, made in natural polypropylene, racked
 DNase, RNase-free, human DNA, PCR inhibitors and Pyrogen free
 STERILE, Autoclavable at 121°C
 Universal filter tips are compatible with most popular brands ,
 Such as: Eppendorf, Gilson, Socorex, HTL, Finnpiptette, Nichiryo, Clinipette brand etc
 All of tips are produced in 100,000net clean room and under strict manufacturing processes
 with premium materials, to assure the highest quality.

Universal Pipette Tips (Filter)



Support Box-packed



Support Refill System



Support Bag-packed

Product Cat#	Description	Color	Volume	Suitable	Packaging Specification
IPT0010-SFRX	Racked, Filter Sterilized	Transparency Blue Yellow	10UL XL	Eppendorf /Gilson/ Thermo	96pcs/Box, 50Box/Carton 4800pcs/Carton
IPT0020-SFR	Racked, Filter Sterilized	Transparency Blue Yellow	20UL	Eppendorf /Gilson/ Thermo	96pcs/Box, 50Box/Carton 4800pcs/Carton
IPT0100-SFR	Racked, Filter Sterilized	Transparency Blue Yellow	100UL	Eppendorf /Gilson/ Thermo	96pcs/Box, 50Box/Carton 4800pcs/Carton
IPT0200-SFR	Racked, Filter Sterilized	Transparency Blue Yellow	200UL	Eppendorf /Gilson/ Thermo	96pcs/Box, 50Box/Carton 4800pcs/Carton
IPT0300-SFR	Racked, Filter Sterilized	Transparency Blue Yellow	300UL	Eppendorf /Gilson/ Thermo	96pcs/Box, 50Box/Carton 4800pcs/Carton
IPT1000-SFR	Racked, Filter Sterilized	Transparency Blue Yellow	1000UL	Eppendorf /Gilson/ Thermo	96pcs/Box, 50Box/Carton 4800pcs/Carton
IPT1250-SFRX	Racked, Filter Sterilized	Transparency Blue Yellow	1250UL / 1000UL XL	Eppendorf /Gilson/ Thermo	96pcs/Box, 50Box/Carton 4800pcs/Carton

Transparen Pipette Tips---Suitable **Beckman**



Product Cat#	Description	Color	Volume	Suitable	Packaging Specification
IPTB20-T-F-S	20ul pipette tips ,filter , sterile ,Transparent box-packed Suitable Beckman	Transparen	20UL	Beckman	96pcs/Box, 50Box/Carton 4800pcs/Carton
IPTB50-T-F-S	50ul pipette tips ,filter , sterile ,Transparent box-packed Suitable Beckman	Transparen	50UL	Beckman	96pcs/Box, 50Box/Carton 4800pcs/Carton
IPTB250-T-F-S	250ul pipette tips ,filter , sterile ,Transparent box-packed Suitable Beckman	Transparen	250UL	Beckman	96pcs/Box, 50Box/Carton 4800pcs/Carton

Transparen Pipette Tips---Suitable Rainin



Universal



Automatic

Product Cat#	Description	Color	Volume	Suitable	Packaging Specification
IPTR20-T-F-S	20ul pipette tips ,filter , sterile ,Transparent box-packed Suitable Rainin Automatic or Universal	Transparen	20UL	Rainin LTS Automatic or Universal	96pcs/Box, 50Box/Carton 4800pcs/Carton
IPTR200-T-F-S	50ul pipette tips ,filter , sterile ,Transparent box-packed Suitable Rainin Automatic or Universal	Transparen	200UL	Rainin LTS Automatic or Universal	96pcs/Box, 50Box/Carton 4800pcs/Carton
IPTR1000-T-F-S	250ul pipette tips ,filter , sterile ,Transparent box-packed Suitable Rainin Automatic or Universal	Transparen	1000UL	Rainin LTS Automatic or Universal	96pcs/Box, 50Box/Carton 4800pcs/Carton

Automatic Pipette Tips---Suitable **Tecan**



Black Conductive



Black Conductive



Clear Non- Conductive

Product Cat#	Description	Color	Volume	Suitable	Packaging Specification
IPTT0050-B-F-S	50ul ,filter , sterile,Black, box-packed Suitable Tecan	Black	50UL	Liquid Handling Arm, Flexible Channel Arm of Tecan platforms	96pcs/Box, 50Box/Carton 4800pcs/Carton
IPTT0200-B-F-S	200ul ,filter , sterile,Black, box-packed Suitable Tecan	Black	200UL	Liquid Handling Arm, Flexible Channel Arm of Tecan platforms	96pcs/Box, 50Box/Carton 4800pcs/Carton
IPTT1000-B-F-S	1000ul ,filter , sterile,Black, box-packed Suitable Tecan	Black	1000UL	Liquid Handling Arm, Flexible Channel Arm of Tecan platforms	96pcs/Box, 50Box/Carton 4800pcs/Carton

Automatic Pipette Tips---Suitable **Hamilton**



50ul Black Conductive



300ul Black Conductive



1000ul Black Conductive

Product Cat#	Description	Color	Volume	Suitable	Packaging Specification
IPTH0050-B-F-S	50ul ,filter , sterile,Black, box-packed Suitable Hamilton	Black	50UL	Hamilton	96pcs/Box, 50Box/Carton 4800pcs/Carton
IPTH0300-B-F-S	300ul ,filter , sterile,Black, box-packed Suitable Hamilton	Black	300UL	Hamilton	96pcs/Box, 50Box/Carton 4800pcs/Carton
IPTH1000-B-F-S	1000ul ,filter , sterile,Black, box-packed Suitable Hamilton	Black	1000UL	Hamilton	96pcs/Box, 50Box/Carton 4800pcs/Carton



**Address – IGNEUS LIFE SCIENCE, 229, TOWER 2, iThink By LODHA
Opp. Xperia Mall, Dombivali East, Thane, 421204**

E-mail: sales@igneus.in Mob: +91 9619000992

Webpage- www.igneus.in