

# U-220

Urine Analyser



## Reliable Testing, Intelligent Design

- User-friendly LCD display and menu-driven interface
- Bright, cold light source for clear results
- Long-lasting performance
- Highly sensitive light sensor with interference resistance- Intelligent pH ratio correction
- Fast, quiet thermal printing
- Seamless computer connectivity via RS-232- Storage for 10,000 test results
- Customizable alert thresholds- Adjustable test sensitivity.
- Easily connect external printers for more printing options.

# U-220

## Urine Analyser



### SPECIFICATIONS

<b>Test Items:</b>	Urobilinogen, Bilirubin, ketone, Creatinine, Blood, Protein, Micro Albumin, Nitrite, Leucocytes, G1ucose, Specific Gravity, pH, Ascorbic Acid, Calcium
<b>Test Color Wave</b>	470nm, 540nm, 620nm
<b>Test Principle</b>	Reflectance and color analysis
<b>Strips Available</b>	10, 11, 14 parameters
<b>Test Speed</b>	120 Tests, 60 Tests/hour and Manually triggered
<b>Display</b>	192 x 64 Graphic LCD display
<b>Data Base Memo</b>	10000 test reports
<b>Printing</b>	Internal 57nm thermal printer
<b>Interface</b>	Standard RS232 interface
<b>USB</b>	Can be connected to the scanner;
<b>Working Environment</b>	Temperature 15-30°C, humidity < 70%
<b>Power Supply</b>	AC220V/110V±15%, 50/60HZ
<b>Power consumption</b>	35VA Size: 253x237x171
<b>Weight</b>	3kg

