



WHY WAIT for test results?



FRENDS™ System



tPSA

**Rapid quantitative assay
for Prostate Specific Antigen**

- Improve patient care
- Increase patient satisfaction
- Increase practice efficiency
- Increase practice profitability

FREND™ System

Accurate, automated testing procedure – simple 3-step operation
No need to batch your sample – run when you are ready

- Controls: Minimum once a month two level liquid controls
- Electronic code chip calibration – no need to run liquid
- Daily electronic system QC control – takes only 1.5 minutes
- Stores 1,000 Test results



Test procedure



STEP 01.
Load sample

>>>



STEP 02.
Insert cartridge

>>>



STEP 03.
Get results


FREND™ PSA Plus

Specification

Assay method	Flourescent immunoassay
Analytical measurement range	0.05 ~ 25 ng/mL
Sample type	Plasma or Serum
Sample volume	35 µL
Limit of detection	0.05 ng/mL
Reaction time	Max 4 min.
1 Box	25 tests
Storage condition	2 °C ~ 8 °C (35 °F ~ 46 °F)

* The specifications are subject to change without prior notice.

Correlation graph



The R^2 correlation value of FREND™ PSA plus is more than 0.96 against laboratory results.

Ordering information

FREND™ System	Cat# F10
FREND™ PSA Plus	Cat# FRPS 025

