

# STW6810-RFID

Sterile Tube Welder



## Generation II

## STW6810-RFID-Sterile Tube Welder

CE MDR Certified / FDA registration NO: BK220673

### Device Specifications

Item Name	STW6810-RFID
Screen Size	10 inch
Dimension	L*W*H:442mm* 760mm* 290mm
Weight	14.5+0.5kg
Working mode	Automatic
Consumable cassette	300 welds per cassette
Voltage	110V-240V

- ✓ Stable Welding Performance
- ✓ Patented Innovation
- ✓ Adjustable Parameters



For a list of worldwide office locations and contact information reach us by E-mail: [info@bms-bloodcare.com](mailto:info@bms-bloodcare.com)

Add : 5 th Floor,C/D,B9,Hi-tech Medical Devices Park,818 Gaoxin Avenue, Wuhan East Lake Hi-tech Development Zone,Wuhan City,Hubei Province

Website : [www.bms-bloodcare.com](http://www.bms-bloodcare.com)

E-mail : [info@bms-bloodcare.com](mailto:info@bms-bloodcare.com)

# STW6810-RFID

Sterile Tube Welder

## 2 Technical Advantages

1 Intelligent welding consumables

2 Digitization of welding information



### Intelligent welding consumables<sup>1</sup>

Built-in RFID chip reader reads the chip on the cassettes in real time, achieving intelligent management on consumables.



### Digitization of welding information<sup>2</sup>

Enables the collection, storage, analysis and export of tube welding data, operator barcode, welding time, blood bag information, weight information etc



**Sterile Tube Welder STW6810-RFID has been exported to 35 countries and regions since 2019, covering Europe, Asia, America, highly accepted by international customers.**



**The first product that integrates welding consumables intelligence and welding information digitization**