

DIRUI

H Urine Analyzer

H-100

www.dirui.com.cn 



Marketed by:

RAPHA [®]
DIAGNOSTICS
Your Trusted Partner In Diagnostics



PRODUCT PROFILE

H-100 Urine Analyzer

Product Characteristics:

- Adopting the advanced high luminosity cold light source with 4 -wavelength, which improves the sensitivity, accuracy, specificity, and reduces the interference from ambient light
- Adopting automatic waste handling equipment, which avoids cross-contamination between samples
- Automatically rectify the test results influenced by non-specificity, pH, specific gravity, and color
- Built-in thermal printer with high speed and low noise; External stylus printer
- Connectable with urine sediment analyzer
- Users can set an abnormal value flag by themselves
- With Chinese and English options
- International, Regular and Symbol system units available for option



Technical Specification:

- Test items: urobilinogen, bilirubin, ketone, blood, protein, nitrite, leukocytes, glucose, specific gravity, pH and VC (or microalbumin)
- Test wavelength: 525nm, 572nm, 610nm, 660nm
- Test principle: Photoelectric colorimetry
- Suitable strips: DIRUI H8, H10 and H11urinalysis strips
- Test throughput:120 strips/h or 60 strips/h optional
- Data memory: 1000 patient results
- Computer interface: RS-232 port; parallel printer interface
- Display: 240×64 dot-matrix LCD
- Language: Chinese and English
- Power supply: 100~240VAC, 50Hz/60Hz
- Power: 40VA
- Dimensions: 376mm×316mm×170mm
- Weight: 3.6kg
- Printer: Built-in thermal printer



Certified to
ISO 9001:2000 and ISO 13485:2003

RAPHA DIAGNOSTICS PVT. LTD.

Corporate Office: 2nd Floor, B-Wing, Vora Plaza, Navghar, Vasai (E), Dist. Thane, Maharashtra - 401 210. INDIA.

Tel.: 91-250-2392924 / 2391740 / 6584304, Helpline No.: 98904 30366, Fax : 91-250-2392878

E-mail : raphadiag@vsnl.net, raphadiag@yahoo.com, Website : www.raphadiag.com